

TdT Assay

pool	per tube		final concentration
25 λ	2.5 λ	1M K ⁺ cacodylate, pH 7.0	100mM
5 λ	0.5 λ	50mM CoCl ₂	1mM
10 λ	1 λ	1:100 dilution of 0.5M dTT	0.2mM dTT
70 λ	7 λ	poly d(A) ₁₀₀ 50	8 8 μ M 3(OT) .0144
50 λ	5 λ	5mM dGTP (1:10 in label)	1mM dGTP
	9 λ	H ₂ O (3H)	
	1 λ	enzyme, straight or diluted	

should be 50mM

TdT#1 (many very small vials)TdT#2 (large screw top vial)amt
enzyme
added

.01 λ \Rightarrow "1-.01"
 .1 λ \Rightarrow "1-.1"
 1 λ \Rightarrow "1-1"

.01 λ \Rightarrow "2-.01"
 .1 λ \Rightarrow "2-.1"
 1 λ \Rightarrow "2-1"

1. Reaction incubated at 37°C for one hour

SR CH 12 1 TIMES
 CR 1 TIMES
 SCR =OFF
 ABC =OFF
 SF =ON
 CALC= 1
 PSI = 10.00 MIN
 CH 1 2.00 2 SIGMA %
 .0 LSR
 .0 BKG
 .00 2 SIGMA B
 0 LL
 397 UL
 CH 2 2.00 2 SIGMA %
 .0 LSR
 .0 BKG
 .00 2 SIGMA B
 655 LL
 940 UL

POS	CH 1	PSZ	CH 2	PSZ	TIME
150 1-1	823.2	2.2	10.3	19.6	10.00
151 1-1	748.7	2.3	9.6	20.3	10.00
152 1-01	710.8	2.3	11.2	18.8	10.00
153 2-1	542.2	2.7	13.7	17.0	10.00
154 2-1	424.0	2.0	11.5	18.6	10.00
	623.4	2.5	21.8	13.5	10.00

MH 68

These calculations have many errors, see p. 10

004

$$\rightarrow \frac{2000 \text{ } \mu\text{g/ml}}{1 \text{ } \mu\text{g/ml}} = 2000 \text{ } \mu\text{g/l}$$

true molarity of DNA

$$A(\text{mw}) \approx 342 \text{ g/mole}$$

$$\left(\frac{2.5 \mu\text{g} \times 10^6 \text{ g}}{\text{mol}} \right) \left(\frac{1}{8 \times 10^3 \text{ mol}} \right) \quad \text{poly(A)}_{100} \text{ (mw)} = (342)(100)$$

$$\frac{(100)(342 \text{ g/mole})}{3.60 \times 10^3 \text{ mole}} \quad \text{should be 50}$$

$$\frac{3.60 \times 10^{-13} \text{ moles}}{25 \times 10^6 \text{ l}} = 1.44 \times 10^{-8} \text{ M} = .044 \mu\text{M}$$

G's added

$$\left(1 \times 10^{-3} \frac{\text{M}}{\text{l}} \right) (25 \times 10^6) N = \# \text{ molecules of G} \quad N = \text{Avogadro's number}$$

$$\left(\frac{1.44 \times 10^{-8} \text{ M}}{1.44 \times 10^{-8}} \right) (25 \times 10^6) N = \# \text{ molecules of 3' OH ends}$$

$$9 \times 10^{-3}$$

$$\text{max possible incorporation of G} = \frac{(1 \times 10^{-3})(25 \times 10^6) N}{(1.44 \times 10^{-8})(25 \times 10^6) N} = 69,444 \text{ G's per 3' end}$$

data here uses 1-1-1 and 1-1-2

$$\frac{800 \text{ cpm}}{(1360 \text{ cpm})(50)} = .0118$$

↑ ↑
pool dot dil. of pool spot

$$(.0118)(69,444) = 819$$

if 2000 cpm

$$\frac{2000 \text{ cpm}}{(1360)(50)} = .0294$$

$$(.0294)(69,444) = 2042$$

$$\left(\frac{1}{50} \right) ($$

0021 T=010.00 A=000094.0(1.5%) B=000255.8(1.5%) C=000200.7(1.5%) R=1.2
 0022 T=010.00 A=000077.2(3.0%) B=000072.7(3.0%) C=0001030.8(2.0%) R=1.2
 0023 T=010.00 A=000151.4(7.0%) B=000187.2(5.0%) C=000214.1(5.0%) R=1.2
 0024 T=010.00 A=000048.8(10%) B=000060.9(10%) C=000083.6(7.0%) R=1.2
 0025 T=010.00 A=000034.4(15%) B=000042.6(10%) C=000069.8(10%) R=1.2
 0026 T=010.00 A=000029.5(15%) B=000039.2(15%) C=000063.7(10%) R=1.34
 0027 T=010.00 A=000028.1(15%) B=000033.0(15%) C=000059.1(10%) R=1.17
 0028 T=010.00 A=000031.9(15%) B=000041.3(10%) C=000072.2(10%) R=1.32
 0029 T=010.00 A=000027.7(5.0%) B=000042.5(3.0%) C=000070.7(3.0%) R=1.25
 0030 T=010.00 A=000055.3(5.0%) B=000039.5(5.0%) C=000062.1(5.0%) R=1.20
 0031 T=010.00 A=000023.2(15%) B=000040.7(15%) C=000056.2(10%) R=1.30
 0032 T=010.00 A=000025.7(15%) B=000037.0(15%) C=000059.3(10%) R=1.36
 0033 T=010.00 A=000030.1(15%) B=000026.8(15%) C=000061.7(10%) R=1.20
 0034 T=010.00 A=000026.1(15%) B=000034.2(15%) C=000057.2(10%) R=1.39
 0035 T=010.00 A=000028.8(15%) B=000036.4(15%) C=000061.7(10%) R=1.38
 0036 T=010.00 A=000030.9(15%) B=000039.6(15%) C=000065.8(10%) R=1.40
 0037 T=010.00 A=000112.5(7.0%) B=000140.9(7.0%) C=000167.7(5.0%) R=1.2

0038 T=010.00 A=000056.6(7.0%) B=000114.9(7.0%) C=000133.7(7.0%) R=1.2
 0039 T=010.00 A=000025.0(15%) B=000036.6(15%) C=000055.7(10%) R=1.2
 0040 T=010.00 A=000022.3(15%) B=000031.6(15%) C=000054.1(10%) R=1.2
 0041 T=010.00 A=000024.4(15%) B=000035.6(15%) C=000058.8(10%) R=1.45
 0042 T=010.00 A=000030.3(15%) B=000042.9(10%) C=000065.5(10%) R=1.40
 0043 T=010.00 A=000020.5(15%) B=000031.7(15%) C=000056.2(10%) R=1.50
 0044 T=010.00 A=000022.7(15%) B=000033.0(15%) C=000056.4(10%) R=1.50
 0045 T=010.00 A=000025.4(15%) B=000034.7(15%) C=000059.1(10%) R=1.35
 0046 T=010.00 A=000026.5(15%) B=000039.8(15%) C=000064.4(10%) R=1.50
 0047 T=010.00 A=001352.4(2.0%) B=001717.4(2.0%) C=001738.5(2.0%) R=1.2
 0048 T=010.00 A=001362.6(2.0%) B=001712.0(2.0%) C=001733.7(2.0%) R=1.2

Corrected for error on p. 2

Unit Calculation - Lang's Way

$$1 \text{ mM dGTP} = \frac{10}{1000} \text{ M} / \lambda$$

10 λ of a 1:100 dilution of pool = 0.1 λ spotted

this produced 1360 counts

$$\frac{13600 \text{ cpm/nucleotide}}{1,360} = 10 \text{ cpm/nucleotide}$$

Assayed highest (4+2) $\rightarrow \frac{2000 \text{ cpm}}{13,600 \text{ cpm/nucleotide}} = \frac{1.5}{1,360} \text{ nmole incorporated}$

$$\frac{1.5}{1,360} \text{ units} / \lambda$$

\rightarrow

$$\frac{1360 \text{ counts}}{0.1 \lambda} \cdot \frac{10 \text{ mM}}{\lambda}$$

$$1360 \frac{\text{cpm}}{0.1 \lambda} \cdot \frac{\lambda}{10 \text{ mM}} = 1360 \frac{\text{cpm}}{\text{mM}}$$

Pucl9 (1 μ g/1 λ)

Sph I \equiv P-S

5 λ Pucl9 (= 5 μ g)
2 λ Sph I (= 10 units)
25 λ 10X low salt buffer
15.5 λ H₂O
25 λ total

Hind III \equiv P-H

5 λ Pucl9 (= 5 μ g)
2 λ Hind III (= 8 units)
25 λ 10X "core" buffer
15.5 λ H₂O
25 λ total

\equiv R-S, .6 μ g/1 λ

1 λ λ DNA (600 μ g/ml)
1 λ Sph I (= 5 units)
2.5 λ 10X L. salt buffer
20.5 λ H₂O
25 λ total

\equiv R-H


1 λ λ DNA
1 λ Hind III (= 4 units)
2.5 λ 10X "core" buffer
20.5 λ H₂O
25 λ total

37° at 10:15

gel

TBE (diluted from 10X buffer)


.7% agarose




P-H R-S R-H

R-S

2/9/84 ~~Atkinson~~ PUC-9H PUC-Hind III



2/10/84 9H1 H12 Puc Alone



Phenol extraction & EtOH ppt of DNA

1. dilute DNA to 200 λ w/ TE
2. Equilibrate Phenol w/ TE & adjust pH to 8.0 w/ 10N N
3. Add 200 λ (equal volume) phenol
4. Spin & suck off phenol (bottom of 2 phases)
5. Spin again & again suck off phenol
6. add 500 λ diethyl ether and let evaporate in hood (but mix first) to remove leftover phenol
7. ppt w/ 22 λ 3M NH_4 acetate
530 λ of 95% EtOH
8. ppt @ -70°C , pellet 5' in mfg, drain & wash
w/ cold 70% EtOH and vac dry

eluted in 15 λ TE

TdT Assay using Puc 19

use	pool	per tube		
100	250 λ	2.5 λ	cacodylate	Pool = 97 λ /tube
	50 λ	0.5 λ	CoCl ₂	
	100 λ	1 λ	1:100 dilution of Ditt	
	500 λ	5 λ	dGTP (1:10 in lab)	
		1 λ	enzyme, undiluted	
450 dGTP		27 λ	DNA (= .8 μ g) @ .27 μ g 1 λ	this amt is roughly scaled up from the puc 100 (multiplied by 500 = 2 x 25) since this DNA is vds \approx about 25x heavier - but it does have 2 ends upon which to add
500 dGTP		13 λ	H ₂ O	
		25 λ total		

2 trials each of puc/H₂O and puc/SPH1
 + control of pool + H₂O

PROG 6 USER 1

SR CH 12 1 TIMES
CR 1 TIMES

SCR =OFF

AGC =OFF

SF =ON

CHLC= 1

PST = 10.00 MIN

CH 1 2.00 2 SIGMA Z

.0 LSR

.0 BKG

.00 2 SIGMA B

0 LL

397 UL

CH 2 2.00 2 SIGMA Z

.0 LSR

.0 BKG

.00 2 SIGMA B

655 LL

940 UL

POS	CH 1	2S2	CH 2	2S2	TIME
-----	------	-----	------	-----	------

139 HindIII-1	27.7	12.0	9.2	20.8	10.00
---------------	------	------	-----	------	-------

140 HindIII-2	29.0	11.7	9.2	20.8	10.00
---------------	------	------	-----	------	-------

141 SphI-1	32.8	11.0	7.5	23.0	10.00
------------	------	------	-----	------	-------

142 SphI-2	35.2	10.6	9.5	20.5	10.00
------------	------	------	-----	------	-------

143 Control	21.5	13.6	11.4	18.7	10.00
-------------	------	------	------	------	-------

144 Pool	4575.3	.9	7.7	22.7	10.00
----------	--------	----	-----	------	-------

→ 107g 1:100 dilution

TdT Assay

97 Pool

77 pcy(A) 100

17 enzyme

87 H₂O

25X

as positive
control

97 Pool

17 enzyme HindIII

47 puc

117 H₂O SphI

25X

→ 10:30pm

Run stopped after 1/2 hour and ppt O/N in
coldroom

TdT Assay

<u>pool</u>	<u>per tube</u>		<u>final concentration</u>
100 λ	25 λ	cacodylate	100mM
20 λ	0.5 λ	CoCl ₂	1mM
40 λ	1 λ	1:100 dilution of Ditt	0.2mM
280 λ	7 λ	poly(A)	8mM 3'OH, 0.4%
200 λ	5 λ	dGTP (1:10 in label)	1mM dGTP
	8 λ	H ₂ O	
	1 λ	enzyme, straight & serially diluted	

3 batches of enzyme were tried at 1 λ , .1 λ , .01 λ , and .001 λ each. Batch 1 was already sorted into small vials. Batch 2 was in the large tube w/ a cap. Vial 3 was uncapped and covered w/ parafilm.

"A-B-C" \Rightarrow A = batch # (1, 2 or 3)

B = quantity of enzyme used (1 λ , .1 λ , .01 λ , or .001 λ)

C = trial number (1 or 2, two trials of each possible used)

C1 = control #1 } pool w/o enzyme in same amt
C2 = control #2 } per tube (25 λ - 1 λ = 24 λ)

pool 1 } 10 λ of a 1:100 dilution of the pool spotted
pool 2 } on a filter and dried

TdT Assay

~~poly(A)~~
 97 pool
 77 poly(A) 100
 17 ~~enzyme~~
~~17~~
~~17~~

Puc/HII

97 pool
 27 DNA @ .27 μ g/17
 57 enzyme
~~17~~
 167

Puc/SphI

97 pool
 27 DNA
 57 enzyme
 167

 λ HII

97 pool
 27 DNA
 57 enzyme
 167

not purified
 in range
 expected

270 pm \rightarrow 50' before stopping

A=000150.6 A=000188.6 (<15%) B=000144.6 (<15%) C=000166.6 (<15%) S=1.050

A=000150.6 A=000184.0 (>20%) B=000192.0 (>20%) C=000222.0 (>20%) S=1.000


A=000150.6 A=000152.0 (>20%) B=000153.0 (>20%) C=000202.0 (<20%) S=7.275

A=000150.6 A=000065.0 (>20%) B=000075.0 (>20%) C=000098.0 (>20%) S=0.975

Control

The above rxn used poly(A)₁₀₀ as a positive control (which worked \Rightarrow the problem is in the DNA and not the ~~the~~ terminal transferase)

Goal to find out if DNA present



N/A
Puc

"Puc"
Puc/HIII

Puc/sphI

\rightarrow from
unreliable
source

\Rightarrow DNA
present

TdT Assay

50x as much				final concentration
pool	per tube			
125 λ	2.5 λ	1 M KTCacodylate pH 7.0	100 mM K ⁺ Ac	
25 λ	0.5 λ	50 mM CoCl ₂	1 mM CoCl ₂	
50 λ	1 λ	1100 dilution of 0.5 M dIT	0.2 mM dIT	
25 λ	0.5 λ	5 mM dGTP	1 mM dGTP	
250 λ	5.0 λ	.087 mM γ -dGTP	.017 mM γ -dGTP	
25 λ	1 λ	.8 μ g/1 λ DNA (= 8 μ g)	1.64 $\times 10^{-7}$ M	
25 λ	1 λ	enzyme		
25 λ	13.5 λ	H ₂ O		
			dGTP 3' ends	= 6,000

HindIII digest of Puc13

90 λ Puc13 DNA (1 μ g/1 λ \Rightarrow 90 μ g)
 22 λ HindIII (= 88 units)
 12.4 10x Cae buffer
~~24.4~~
 124.4 λ total

4:30 pm


Digest was checked after 5 hrs, but then allowed to go on

30 Minutes

David's Engine (Sec 28 units A)

poly (CH₂)₁₀ 1.000.50 A=0001 68.00 (20%) B=0001 55.00 (20%) C=0000 202.00 (20%) R=1.000
poly (CH₂)₁₀ 1.000.50 A=0001 68.00 (20%) B=0001 55.00 (20%) C=0001 202.00 (20%) S=1.000

after 4 hours



12/epH
 Puc19 nucA
 Puc13 nucA
 12/HindIII

TdT Assay

<u>Puc13</u>	<u>Puc13/HindIII</u>	<u>Puc13/HindIII + Puc19</u>	<u>Control</u>
97 pool	97 pool	97 pool	97 pool
77 poly(A)	1000 DNA	1000 DNA	
17 enzyme	17 enzyme	77 poly(A)	
87 H ₂ O	14.57 H ₂ O	17 enzyme	
257	257	7.57 H ₂ O	
		257	

→ 30 min.

100.50 A=002976.0(7.0%) B=003714.0(5.0%) C=003770.0(5.0%) R=1.250
 100.50 A=000118.0(>20%) B=000123.0(>20%) C=000158.0(>20%) R=1.400
 100.50 A=004618.0(5.0%) B=005926.0(5.0%) C=005984.0(5.0%) R=1.236

#G's per ~~1000~~ 3' end of Puc13

Molarity of Puc13:

$$\frac{(1 \times 10^{-6} \text{ g}) (1 \text{ g})}{1.8 \times 10^6 \text{ g/mole}} = 5.56 \times 10^{-13} \text{ moles (in 120 } \mu\text{L DNA)}$$

$$\frac{5.56 \times 10^{-13} \text{ moles}}{0.025 \text{ mL}} = 2.22 \times 10^{-11} \text{ M}$$

there are 2 3' ends per molecule $\Rightarrow (2)(2.22 \times 10^{-11} \text{ M})$
 $= 4.44 \times 10^{-11} \text{ M}$

(cf. Wu who recommends
 $\frac{60 \text{ pmole}}{\text{mL}} = 6 \times 10^{-11} \text{ M}$)

$$\frac{5.56 \times 10^{-13} \text{ moles}}{25 \times 10^{-6} \text{ L}} = 2.22 \times 10^{-8} \text{ M} = 2.22 \times 10^{-5} \text{ mM}$$

molarity of GTP

remember
 2 3' ends per molecule $\Rightarrow 4.44 \times 10^{-5} \text{ M (3'}$

$4.5 \Rightarrow 25 \Rightarrow \frac{4.5}{25} = .18 \text{ dilution}$

$(.18)(50 \text{ mM}) = 9 \text{ mM}$

max possible
incorp

$$\frac{GTP}{PuCl_3} = \frac{9 \mu M}{2.22 \times 10^{-5} M} = 405,405 \text{ G's/PuCl}_3$$

$$\frac{GTP}{PuCl_3} = \frac{9 \mu M}{4.44 \times 10^{-5} M} = 202,703 \text{ G's/3'OH}$$

fraction of counts incorporated
background

$$\frac{118 - 30}{(250 \times 4575) \text{ cpm}} = \frac{88}{(250 \times 4575)} = 7.7 \times 10^5$$

.17 g
spooled
pool

4575 cpm

in 257 g pool
there are
(250 x 4575 cpm)

actual
incorporation:

$$7.7 \times 10^5 (202,703) = 156 \text{ G's/3' end}$$

#G's per 3' end on poly(A)

molarity of poly(A) (see also p. 4)

25 g

$$\frac{(2.5 \times 10^{-6} \text{ g}) (7 \times 10^{-3} \text{ mol})}{(50 \times 347 \text{ g/mol})} = 5.04 \times 10^{-11} \text{ mole}$$

$$\frac{5.04 \times 10^{-11} \text{ mole}}{25 \times 10^{-6} \text{ l}} = 2.02 \times 10^{-6} M$$

$$= 2.02 \times 10^{-3} \text{ mM}$$

max possible incorp.

$$\frac{GTP}{PuCl_3} = \frac{9 \mu M}{2.02 \times 10^{-3} M} = 4,455 \text{ G's/3'OH}$$

Notice that although $\frac{G}{3'end}$ is about the same for both Puc13/HIII and poly(A)₅₀ the ratio of $\frac{GTP}{3'OH}$ is greater for Puc13 (\Rightarrow equilibrium-mass action-effects should be stronger)

Since $\frac{GTP}{3'OH}$ for Puc13 is about 50X greater than the ratio for poly(A)₅₀ and the amt of incorp is about the same, this suggests the Puc13 rxn is 50X less efficient than poly(A)₅₀

(Obvious: The reason why there is relatively few counts depicted for Puc13/HIII is because ~~the~~ the molarity of this is much less than that of the poly(A)₅₀ as used in the same assay. To better differentiate background from incorporation in the Puc rxn, I should ↑ the amt of label in the pool-by sens, 10X (sounds good).)

Fraction actually incorporated

$$\frac{30000}{(250)(4575)} = .0026$$

actual
incorp $(.0026 \times 4,455) = \boxed{12 \text{ (6/3) ed}}$

The molarity of $^3\text{H dGTP}$ in stock

plate
from
Spec →
sheet

$$\frac{.013 \times 10^{-3} \text{ g}}{(.25 \times 10^{-3} \text{ L})(595.2 \text{ g/mol})} = 8.7 \times 10^{-5} \text{ M} = 87 \times 10^{-6} \text{ M} = .087 \text{ mM}$$

poly(A)so
9λ pool
7λ poly(A)
1λ enzyme
8λ H₂O
25λ total

Luc13/H^{III}
9λ pool
1λ DNA
1λ enzyme
14λ H₂O
25λ total

Control
9λ pool
2λ DNA (poly(A)so)
1λ enzyme
8λ H₂O
20λ total

put in at 37° at 3:30pm
→ allowed to go for 18 hours

Poly(A)so

015 T=000.50 A=000413.0(5.0%) B=005312.0(5.0%) C=005369.0(5.0%)
016 T=000.50 A=000178.0(>20%) B=000176.0(>20%) C=000300.0(>20%)
017 T=000.50 A=000052.0(>20%) B=000070.0(>20%) C=000104.0(>20%)
18 WASH ASSAY

$\text{PVC}13/\text{H}^{\text{III}} \Rightarrow 226/3^{\text{rd}}$

$\text{payd}(\text{A})50 \Rightarrow 176/3^{\text{rd}}$

~~dGTP, Na^+ salt~~

~~Pl. in new has .025g~~

~~507.20~~

~~609 $\frac{\text{g}}{\text{mole}}$~~

~~1ml 1M solution
 \hookrightarrow~~

TdT Assay

This assay uses $\frac{1}{10}$ the molar concentration of Pvc13/HIII and $\frac{1}{100}$ as the previous series. This effectively \uparrow the relative dGTP/3'OH ratio by a factor of 10. Additionally the label has been increased by 14 fold. (There should be more incorporation per 3' end, but there are fewer 3' ends, so this increase is not expected to appreciably alter the counts, relative to background.)

1. $\frac{1}{10}$ the previous concentration of DNA $\Rightarrow 2.03 \times 10^6$ excess of dGTP over 3'

✓ 17 of a 1:10 dilution of μ g17 Pvc13/HIII DNA

$\cdot 1 \mu\text{g}$
 $= 4.44 \times 10^6 \text{ mM}$

17 TdT

3 units

✓ 57 50 mM dGTP

100 mM

147 ^3H dGTP

$\cdot 05 \text{ mM}$

✓ Pool \equiv { 2.57 1M K⁺ cacodylate

100 mM

547 { 0.57 50 mM CoCl₂

1 mM

17 1:100 dilution of 0.5M dIT

0.2 mM dIT

257 total

2. ~~17 of a 1:100 dilution of~~

$\frac{1}{100}$ the previous concentration of DNA $\Rightarrow 2.03 \times 10^7$ excess of dGTP over 3'

17 of a 1:100 dilution of μ g17 Pvc13/HIII DNA

$\cdot 01 \mu\text{g}$
 $= 4.44 \times 10^{-7} \text{ mM}$

Because the rest is the same as the above assay

257

3. As before,

17 of μ g17 Pvc13/HIII DNA

The rest is the same as above

4. Control, everything, but no DNA or enzyme

2.03×10^5 $\frac{\text{cpm}}{3' \text{ end}}$

T₁ 17 of the 237 total control rxn was spotted for use as a "part spot"

P₁₀₀ ⊖ ← signifies no dGTP
(enough for 100 257 assays)

2507 cacodylate

507 CoCl₂

997 H₂O

17 undiluted 0.5M dIT

2 parts were prepared - red label CoCl₂ added before dIT (added last) - turned brown - not used. When CoCl₂ was added last there were no brown ppt. This good stuff is in a green labelled tube and was used for the above rxns.

11-23-78 T=000.50 A=000670.0(15%) B=000732.0(15%) C=000754.0(15%) R=1.090
11-23-78 T=000.50 A=000334.0(20%) B=000456.0(15%) C=000424.0(15%) R=1.157
11-23-78 T=000.50 A=000400.0(15%) B=000430.0(15%) C=000512.0(15%) R=1.200
11-23-78 T=000.50 A=000460.0(15%) B=000558.0(15%) C=000570.0(15%) R=1.173
11-23-78 T=000.50 A=211804.0(0.7%) B=248908.0(0.7%) C=250060.0(0.7%) S=1.174

Perhaps the excess etOH inhibited the rxn, or there may also be substrate level inhibition

of 1000000000
 1000000000
 1000000000
 1000000000

the other side is the fact that the system is not
 - related to the fact that the system is not
 that the system is not related to the fact that the
 system is not related to the fact that the system is not
 related to the fact that the system is not related to the
 fact that the system is not related to the fact that the
 system is not related to the fact that the system is not

★
Idea: put in *dTPP as trace
 in dGTP failure. What will
 this tell?

30 minutes

TdT assay

1. 17 μ g DNA-Puc¹³H^{III} 2 μ g (17) 1 μ g
 17 TdT 3 units
 17 α -³²P dTTP (Evaporated and
 resubstituted)
 47 Pool- 51d
 187 50 mM dTTP 36 mM dTTP
 257

2. Same as above
 but 17 μ g 1 μ g (17) Puc¹³H^{III} 1 μ g

3. as above, but 17 μ g 0.1 μ g (17) Puc¹³H^{III} 0.1 μ g

4. control
 no DNA
 no TdT

237 total

~~spot on gel~~

dilute 17 1:10 then spot

- 5-5. Check to see that such massive units of
 dNTP's are all causing substrate-level inhibition

17 μ g 1 μ g (17) Puc¹³H^{III}

17 TdT

17 α -³²P dTTP

47 Pool-

17 ~~Pool~~ of 1:4 dilution of 50 mM dTTP \leftarrow 5 mM

g. a 1:4 dilution of 50 mM dTTP

177 H₂O
 257 totalPool Spot: 17 μ g 1:10 dilution
 from control (237 total)

PROG 6 USER 1

SR CH 12 1 TIMES

CR 1 TIMES

SCR =OFF

RQC =OFF

SF =ON

CALC= 1

PST = 10.00 MIN

CH 1 2.00 2 SIGMA Z

.0 LSR

.0 BK6

.00 2 SIGMA B

0 LL

397 UL

CH 2 2.00 2 SIGMA Z

.0 LSR

.0 BK6

.00 2 SIGMA B

655 LL

940 UL

	POS	CH 1	2SZ	CH 2	2SZ	TIME
1	8	156.4	5.0	5622.1	.8	10.00
10	9	22.9	13.2	472.8	2.9	10.00
10	10	22.7	13.2	316.9	3.5	10.00
S	11	3624.6	1.9	202151.8	.2	2.76
C	12	18.4	14.7	235.9	4.1	10.00
par	13	281.3	3.7	9136.4	.6	10.00

Calculations for above experiment

1. 1700 1 μg/λ DNA

$$\frac{GTP}{PUC13} = \frac{10 \text{ mM}}{4.44 \times 10^5 \text{ mM}} = 2.25 \times 10^5 \leftarrow \text{max possible incorporation}$$

Fraction of counts incorporated

$$\hookrightarrow \frac{5622.1 - 235.9}{(9136.4)(1.087 \times 10^2)} = \frac{5386.2}{2,482,816.7} = .00217$$

↑ diluted 1:10
x 25 to get total possible counts

$$\left(\frac{1}{23}\right)(25\lambda) = 1.087\lambda$$

↑
1λ from a total volume of 23λ is equal to 1.087 from a final volume of 25λ

actual incorporation

$$(.00217)(2.25 \times 10^5) = 488 \text{ T/3' end}$$

2. 1700 1 μg/λ DNA

$$\frac{GTP}{PUC13} = \frac{10 \text{ mM}}{4.44 \times 10^5 \text{ mM}} = 2.25 \times 10^5 \leftarrow \text{max possible incorp}$$

Fraction of counts incorporated

$$\frac{478.2 - 235.9}{2,482,816.7} = 9.76 \times 10^{-5}$$

actual incorporation $\rightarrow (9.76 \times 10^{-5})(2.25 \times 10^5) = 220 \text{ T/3' end}$

B.17 of .01 $\mu\text{g/l}$ Pu¹³⁴

$$\frac{\text{max possible}}{\text{incorp}} \rightarrow 2.25 \times 10^7$$

Fraction of counts
incorp

$$\rightarrow \frac{3169 - 235.9}{2,482,816.7} = 3.26 \times 10^{-5}$$

$$\frac{\text{actual}}{\text{incorp}} \rightarrow (2.25 \times 10^7) (3.26 \times 10^{-5}) = 734$$

~~Q28~~

S.E.S 1700 $\mu\text{g/l}$ Pu¹³⁴ - diff conc of dTIP

$$\text{vicinity of dTIP} \left(\frac{1}{100}\right) (50\text{mm}) = .5\text{mm}$$

$$\frac{\text{max possible}}{\text{incorp}} = \frac{\text{dTIP}}{3'01} = \frac{.5\text{mm}}{4.44 \times 10^5\text{mm}} = 1.13 \times 10^4$$

Fraction of counts
incorp

$$\rightarrow \frac{202,151.8 - 235.9}{2,482,816.7} = .0813$$

$$\frac{\text{actual}}{\text{incorp}} \rightarrow (.0813) (1.13 \times 10^4) = 919 \text{ T/3' end}$$

TdT Assay

1. 47 part
 17 Puc13/H^{III} @ 1 μ g/17
~~100~~ 57 ³H dGTP - evap & reconstituted
 17 TdT 1:10
 57 ~~1:10~~ dilution of 50 mM dGTP
 97 H₂O
257

std
 1 μ g
 ~3 units
 1 mM

2. 47 part
 17 Puc13/H^{III}
 57 ³H dGTP - evap & reconstituted
 17 TdT
 107 1:10 dilution of 50 mM dGTP 2 mM dGTP
 47 H₂O
257

3. Same as above, but
 2.57 1:10 dil. of 50 mM dGTP
 11.57 H₂O
257

0.5 mM dGTP

4. 17 1:10 dil. of 50 mM dGTP
 137 H₂O
257

0.2 mM dGTP

2000 of ³H dGTP was vac dried & reconstituted
 in 407 TE (\Rightarrow 5X concentration)

✓ 100 μ l H_2O ✓ 10 μ l Puc 13/HIII @ 100 μ g/l100 μ g
~3 units10 μ l TdT10 μ l α - ^{32}P dTTP✓ 40 μ l pool✓ 20 μ l of a 1:4 dilution of 50 mM dTTP \Rightarrow 1 mM dTTP
run for 25 μ l total

10 min

30 min

1 hour

2 hours

4 separate runs

NO DNA ~~over~~ TdT

+ control (230 total)

+ pool-spot 10 from control

+ "Save" 30 min for the sequencing experiment

On "Save" - don't stop w/ acid wash-step
by chilling - it will later be phenol extracted
& then cut & run out on gel

10 min 005 T=000.50 A=258594.0(0.7%) B=297376.0(0.7%) C=303772.0(0.7%) R=1.151

30 min 006 T=000.50 A=273180.0(0.7%) B=275920.0(0.5%) C=332902.0(0.5%) R=1.193

1 hour 007 T=000.50 A=300768.0(0.7%) B=353172.0(0.5%) C=365600.0(0.5%) R=1.190

2 hours 008 T=000.50 A=195122.0(0.7%) B=219522.0(0.7%) C=223372.0(0.7%) R=1.125

009 T=000.50 A=000664.0(1.5%) B=000768.0(1.5%) C=000856.0(1.0%) R=1.151

010 T=000.50 A=151398.0(1.0%) B=200524.0(0.7%) C=207264.0(0.7%) S=1.151

could pool 12 addl. from orig. 230 control

Remember, the calculated values are per 3' end

30 minutes

$$\frac{dTTP}{dT/3'end} = \frac{1mM}{4.44 \times 10^{-5}mM} = 2.25 \times 10^4$$

$$\frac{273,180 - 664}{(151,398) \left(\frac{25}{23}\right) (25)} = \frac{272,516}{4,114,076.1} = .0662$$

$$(.0662)(2.25 \times 10^4) = 1,490 T/3'end$$

10 minutes

$$\frac{258,694 - 664}{4,114,076.1} (2.25 \times 10^4) = 1,411 T/3'end$$

1 hour

$$\left(\frac{300,768 - 664}{4,114,076.1} \right) (2.25 \times 10^4) = 1,641 T/3'end$$

2 hours

$$\left(\frac{195,122 - 664}{4,114,076.1} \right) (2.25 \times 10^4) = 1,064 T/3'end$$

check to see if α - 32 P dTTP

-phenol/extract the sample for the gel
-try run 4 in labelling rxn

idea: in plastic bag put
5 μ mol of all 4 bases

TdT Assay w/ all 4 bases present

1.	1 μ l	Puc 13/H ^{III} DNA @ 1 μ g/l	1 μ g
	1 μ l	TdT @ ~30 U/l	30
	4 μ l	Pool-	std
	1.3 μ l	1:10 dil of 50 mM dGTP (= 5 mM)	0.25 mM dGTP
	1.3 μ l	1:10 dil of 50 mM dCTP (= 5 mM)	0.25 mM dCTP
	1.3 μ l	1:10 dil of 50 mM dATP (= 5 mM)	0.25 mM dATP
	1.3 μ l	1:10 dil of 50 mM dTTP (= 5 mM)	0.25 mM dTTP
	1 μ l	α - 32 P dTTP (evap & reconstituted)	32 P-T-label
	5 μ l	3 H dGTP (evap & reconstituted)	3 H-G-label
	7.8 μ l	H ₂ O	
	25 μ l	total	

2.	1 μ l	Puc 13/H ^{III} DNA @ 1 μ g/l	1 μ g
	1 μ l	TdT @ ~30 U/l	30
	4 μ l	Pool-	std
	1.3 μ l	1:10 dil of 50 mM dGTP	0.25 mM dGTP
	1.3 μ l	" dCTP	" dCTP
	1.3 μ l	" dATP	" dATP
	1.3 μ l	" dTTP	" dTTP
	1 μ l	α - 32 P dCTP	32 P-C-label
	5 μ l	3 H dATP (evap & reconstituted)	3 H-A-label
	7.8 μ l	H ₂ O	
	25 μ l	total	

G₁ = 3. Control - for rxn #1

4 μ l	Pool-
1.3 μ l	1:10 dil of 50 mM dGTP
1.3 μ l	" dCTP
1.3 μ l	" dATP
1.3 μ l	" dTTP
1 μ l	α - 32 P dTTP
5 μ l	3 H dGTP
9.8 μ l	H ₂ O

G₂ = 4. Control - for rxn #2

4 μ l	Pool-
1.3 μ l	dGTP
1.3 μ l	dCTP
1.3 μ l	dATP
1.3 μ l	dTTP
1 μ l	α - 32 P dCTP
5 μ l	3 H dATP
9.8 μ l	H ₂ O

From each of the controls spot ¹⁷ ~~17~~ as the "pool"

~~1, 2, C1, C2, P1, P2~~

1, 2, C1, C2, P1, P2

→ ~~40 min~~ 40 min
~~20 min~~

INCORP. of G

$$\frac{\text{max possible G incorp}}{\text{G incorp}} \rightarrow \frac{dGTP}{3'end} = \frac{0.25mM}{4.44 \times 10^{-5}mM} = 5,631$$

$$\frac{\text{fraction of counts incorp}}{\text{incorp}} \rightarrow \frac{136,890 - 570}{(122,262)(25)} = .0446$$

$$\frac{\text{actual incorp}}{\text{incorp}} \rightarrow (.0446)(5,631) = 251$$

INCORP. A

$$\frac{\text{max possible A incorp}}{\text{A incorp}} \rightarrow \frac{0.25uM}{4.44 \times 10^{-5}uM} = 5,631$$

$$\frac{139,538 - 584}{(131,524)(25)} (5,631) = 238$$

~~incorp~~

incorp of C

$$\frac{201,442 - 926}{(412,726)(25)} (5,631) = 109$$

incorp of T

$$\frac{47,096 - 266}{(113,821)(25)} (5,631) = 93$$

	base/blend	%
A	238	34.4%
C	109	15.8%
G	251	36.3%
T	93	13.5%
	691	100%

$\Rightarrow G > A > C > T$

T=001.00 A=000570.0(10%) B=000266.0(15%) C=000999.0(7.0%) R=0.456
 T=001.00 A=000534.0(10%) B=000926.0(7.0%) C=001322.0(5.0%) R=1.536
 T=001.00 A=136790.0(0.7%) B=047096.0(1.0%) C=226411.0(0.5%) R=0.342
 T=001.00 A=139532.0(0.7%) B=201402.0(0.5%) C=416207.0(0.3%) R=1.422
 T=001.00 A=122262.0(0.7%) B=112821.0(0.7%) C=302295.0(0.5%) R=0.821
 T=001.00 A=131524.0(0.7%) B=412726.0(0.3%) C=705104.0(0.3%) S=0.921

A: $1/2^3H$ B: $3/2^6H$ C: $3H + 14C + 22P$

run over again, exactly the same, w/ the
 exception of #2, which will just be stopped
 two #2's & run out on a gel

~~3.85~~ 40 minutes, as before

002 T=001.00 A=000959.0(2.0%) B=000357.9(1.5%) C=001562.0(7.0%) R=0.368
 005 T=001.00 A=000549.0(1.0%) B=000516.0(1.0%) C=001333.0(7.0%) R=0.942
 006 T=001.00 A=281265.0(0.5%) B=049662.0(1.0%) C=432974.0(0.3%) R=0.947
 007 T=001.00 A=174665.0(0.5%) B=279508.0(0.5%) C=551349.0(0.3%) R=1.400
 008 T=001.00 A=172395.0(0.5%) B=113914.0(0.7%) C=346891.0(0.5%) R=0.642
 009 T=001.00 A=183277.0(0.5%) B=324967.0(0.5%) C=593820.0(0.3%) R=1.773

same channels as above

⇓

$$\boxed{G} \rightarrow \frac{281,265 - 959}{(177,395)(25)} (5,631) = 356$$

$$\boxed{A} \rightarrow \frac{174,665}{\cancel{183,209} - 549} (5,631) = 214$$

$$(183,277)(25)$$

$$\boxed{C} \rightarrow \frac{279,508 - 516}{(324,967)(25)} (5,631) = 193$$

$$\boxed{T} \rightarrow \frac{69,662 - 357}{(113,914)(25)} (5,631) = 137$$

G	356	<u>9%</u> 40%
A	214	24%
C	193	21%
T	137	15%
Total	900	

$\Rightarrow G > A > C > T$

The saved portion from the peptide tailing (p.15 yellow) was frozen after the rxn was completed (30 min). Then, together w/ the 'saved portion' from the experiment immediately above, the two were phenol extracted & EtOH ppt'd

↓ re-eluted
each, ^{12.5}~~12.5~~ H₂O ← 1 μg in rxn mix
17 spl = 5 units
1.5 10x low salt buffer
15.7 total

37°C @ 8:40 am

↳ 2 hours, 20 mins long

↳ both of these samples are loaded on a sequencing gel run by class

- Calculations of Stop codon appearance based on base composition

Term. codons: UAA UGA \Rightarrow TAA TGA
UAG TAG

Using data from the first experiment w/ all 4 n.t.'s present

$$(T)(A)(A) + (T)(A)(G) + (T)(G)(A)$$

Where $T = \text{prob of } T = .13$
 $A = \text{etc.} = .344$
 $C = .158$
 $G = .363$

$$= (T)(A^2) + 2(T)(A)(G)$$

$$= (.135)(.344)^2 + (2)(.135)(.363)(.344)$$

$$= .04969$$

# of residues	fraction terminated	
10	.95	.983
30	.60	.843
50	.22	.60
75	.078	.427
100	.022	.279
150	.006	.182
200	.00047	.078
300	.000037	.033
	.0000002	.006

$$A = .15$$

$$T = .15$$

$$C = .35$$

$$G = .35$$

$$(.15)(.15)^2 + 2(.15)(.15)(.35) =$$

$$A = .15$$

$$T = .15$$

$$C = .40$$

$$G = .30$$

Anticipated N.T. concentrations

$$A: \left(\frac{.34}{.15} \right) (.25 \text{ mM}) =$$

$$\begin{array}{l} \text{desired} \rightarrow \left(\frac{.15}{.34} \right) (.25 \text{ mM}) = .11 \text{ mM} \\ \text{actual} \rightarrow \end{array}$$

$$T: \left(\frac{.15}{.135} \right) (.25 \text{ mM}) = .28 \text{ mM}$$

$$C: \left(\frac{.40}{.158} \right) (.25 \text{ mM}) = .63 \text{ mM}$$

$$G: \left(\frac{.30}{.363} \right) (.25 \text{ mM}) = \frac{.21 \text{ mM}}{.12 \text{ mM}}$$

1.23 mM total [NT]

these values are calculated per
μg of DNA

TdT Assay (Try out previously calculated concentrations & check the new P-L terminal transferase)

1. 17 ~~0.01~~ ~~0.001~~ ~~0.0001~~ μ g PUC13/H^{III} DNA @ 1 μ g/l
 17 lab TdT (~30 units/l)
 47 pool-
 2.27 1:40 dil of 50mM dATP (=1.25mM)
 1.47 1:10 dil of 50mM dTTP (=5mM)
 3.27 1:10 dil of 50mM dCTP (=5mM)
 1.17 1:10 dil of 50mM dGTP (=5mM)
 17 α -³²P dTTP
 37 ³HdGTP
 7.17 H₂O

~~Exp~~ 1 μ g
 30 units
 std concentration
 .11mM dATP
 .28mM dTTP
 .63mM dCTP
 .21mM dGTP

2. ~~0.01~~ same as above except

17 α -³²P dCTP
 37 ³HdATP
 257

3. Control for rxn 1

as above #1 but no DNA or TdT
 \Rightarrow 9.17 H₂O
 257

4. Control for rxn 2

as #2 but no DNA or TdT
 \Rightarrow 9.17 H₂O
 257

P-L TdT

repeat above: ~~use same control;~~
 1 & 2 use just 0.57 TdT
 and compensate w/ 7.67 H₂O
 \Rightarrow 1' & 2' then labelled




Take 27 from each control and spot to
 one at a "pool" spot

~~#45~~ \rightarrow 40 minutes

Stock ATP was diluted, so that there'd be enough for this experiment

124 T=000.50 A=000292.0(20%) B=000284.0(20%) C=000742.0(15%) S=1.000
 125 T=000.50 A=000166.0(>20%) B=000389.0(10%) C=001338.0(10%) S=5.500
 126 T=000.50 A=221554.0(0.7%) B=040312.0(1.5%) C=350908.0(0.5%) S=0.271
 127 T=000.50 A=013949.0(3.0%) B=494349.0(0.5%) C=593345.0(0.5%) S=3.211
 128 T=000.50 A=031372.0(2.0%) B=007889.0(3.0%) C=049802.0(1.5%) S=0.283
 129 T=000.50 A=002472.0(7.0%) B=072920.0(1.5%) C=033729.0(1.0%) S=9.442
 130 T=000.50 A=230414.0(0.7%) B=273562.0(0.7%) C=634450.0(0.5%) S=1.202
 131 T=000.50 A=029454.0(2.0%) B=692174.0(0.5%) C=949400.0(0.2%) S=1.207

Same channels as before

A = 
 B = 
 C = 

Lab TdI

$$[A] \quad \frac{dATP}{3'} = \frac{11 \text{ mM}}{4.44 \times 10^{-5} \text{ mM}} = 2,478$$

$$\frac{13,948 - 166}{29,454 (12.5)} (2,478) = 93$$

$$[G] \quad \frac{dGTP}{3'} = \frac{.21 \text{ mM}}{4.44 \times 10^{-5} \text{ mM}} = 4,730$$

$$\frac{221,554 - 292}{239,414 (12.5)} (4,730) = 363$$

$$[C] \quad \frac{dCTP}{3'} = \frac{.63 \text{ mM}}{4.44 \times 10^{-5} \text{ mM}} = 14,189$$

$$\frac{494,348 - 888}{199,174 (17.5)} (14,189) = 801$$

$$\boxed{T} \quad \frac{dTP}{31} = \frac{.28 \mu M}{4.44 \times 10^{-5} \mu M} = 6,306$$

$$\frac{60,312 - 284}{278,562(12.5)} (6,306) = 109$$

P-L TdT

$$\boxed{A} \quad \frac{2,472 - 166}{29,454(12.5)} (2,478) = 16$$

$$\boxed{G} \quad \frac{3,372 - 292}{230,414(12.5)} (4,730) = 51$$

$$\boxed{C} \quad \frac{22,920 - 888}{699,174(12.5)} (14,189) = 117$$

$$\boxed{T} \quad \frac{8,898 - 284}{278,562(12.5)} (6,306) = 16$$

		<u>P-L</u>		<u>Lab Prep</u>		<u>Calculated goal</u>
	#		%	#	%	
Lab Prep	A	16	8%	93	7%	15
Lab Prep	C	117	59%	801	59%	40
Lab Prep	T	16	4%	109	8%	15
	G	51	26%	363	27%	30
	total	200		1366		

appearance of term. codons for above base ratios

$$(T)(CA)^9 + 2(CT)(A)^6$$

$$(0.08)(.08)^2 + (2)(.08)(.08)(.26)$$

<u>residues</u>	<u>Procn stopped</u>
1	.9962
50	.825
100	.68
150	.56
200	.46
250	.38
300	.32
350	.26
400	.21
500	.15
600	.10
700	.07
1000	.02

← 455 ⇒ 17% still untermint
(early aa
are high
symmetrical)

⇒ look to see how aa distribution ratio is affected

A New calculation of anticipated [NTP] to get
the desired

$$A = .15$$

(chosen on p. 19 white)

$$T = .15$$

$$C = .40$$

$$G = .30$$

A: $\frac{.15}{.08} (.11) = .24 \text{ mM}$ ← though I know that this ^{conc.} ratio before
gave me too much

T: $\frac{.15}{.08} (.28) = .53 \text{ mM}$

C: $\frac{.40}{.159} (.63) = 1.42 \text{ mM}$

G: $\frac{.30}{.27} (.21) = \frac{.23 \text{ mM}}{1.42 \text{ mM}}$

I think that the problem of this sort of calculation is
that one is dependent upon the other so that these
should be "normalized" w/ respect to each other in any
calculation

desired goal:

$$\begin{aligned} A &= .05 \\ C &= .25 \\ G &= .25 \\ T &= .45 \end{aligned}$$

This is one way to
calculate expected
concentrations

& also sure enough to run out on a gel

fraction in product	fraction in mix	rxn fraction product fraction = M
A = .07	.11 mM	.082/.07 = 1.17
C = .59	.63 mM	.47/.59 = .80
T = .08	.28 mM	.209/.08 = 2.61
G = .27	.21 mM	.157/.27 = .581
	<u>1.34 mM</u>	

product
desired ~~rxn~~ fraction

$$\begin{aligned} A &= .05 \\ C &= .25 \\ G &= .25 \\ T &= .45 \end{aligned}$$

rxn fraction

(estim.) = M (desired product fraction)

$$\begin{aligned} &.0684 \\ &.20 \\ &.653 \\ &\underline{.261} \\ &1.18 \end{aligned}$$



normalize w/ respect to total

A
C
G
T

$$\begin{aligned} &.0684/1.18 = .058 \\ &.20/1.18 = .17 \\ &.653/1.18 = .55 \\ &\underline{.261/1.18 = .22} \\ &.998 \end{aligned}$$

1. 17 μ g 131H DNA @ 1 μ g/17

17 TdT

47 pool

1.27 1:40 aile of 50 mM dATP (= 1.25 mM)

3.47 1:40 aile of 50 mM dCTP (= 1.25 mM)

2.87 1:10 aile of 50 mM dGTP (= 5 mM)

1.17 1:10 aile of 50 mM dTTP (= 5 mM)

17 α - 32 P dTTP

37 3 H dGTP

6.57 H_2O

257

1 μ g
~3 units
Std concentrations
~~0.58 mM dNTP~~ ~~0.05~~
1.7 mM dCTP
5.5 mM dGTP
2.2 mM dTTP

2. Same as above but

17 α - 32 P dCTP

37 3 H dATP

Save \equiv 1

For controls pool: before adding TdT to the 3 rxn mixtures draw out 20 and spot to count - replace w/ 20 of H_2O . So, therey'll be a

C1
C1' (\equiv Save)
C2

Another way to calculate out expected concentrations

Via ratio of dNTP/3'

		$[dNTP/B]$		anti-tail $[dNTP]$
A = .07	.11 μM	2,478	$.05 / .07 (2,478) = 1,770$.079 μM
C = .59	.63 μM	14,189	$.25 / .59 (14,189) = 6,012$.27 μM
T = .08	.28 μM	6,306	$.25 / .08 (6,306) = 19,706$.87 μM
G = .27	.21 μM	4,730	$.45 / .27 (4,730) = 7,883$.35 μM

↑
total tail length
1366

total tail length of 1,000 would be better, so multiply
all $[dNTP]$'s by $\frac{1,000}{1366}$ ($\approx .732$)

($[dNTP]$)
A .058 μM
C .20 μM
T .64 μM
G .26 μM
1.16 μM

1. 81 λ PVL131HIII DNA @ 1 μ g 1 λ

1 λ TdT

4 λ pool

1.2 λ 1:40 a11g sermd ATP

1 λ 1:10 a11g sermd dCTP

3.2 λ 1:10 a11g sermd TTP

1.3 λ 1:10 a11g sermd dGTP

1 λ α -32p dTTP

3 λ 3 H dGTP

8.2 λ H_2O

25 λ total

1 μ g
~3 units

0.5 μ g dATP	0.06
0.5 μ g dCTP	0.06
0.5 μ g dTTP	0.06
0.5 μ g dGTP	0.06
0.5 μ g sermd ATP	0.06
0.5 μ g sermd dCTP	0.06
0.5 μ g sermd TTP	0.06
0.5 μ g sermd dGTP	0.06

2. Same as above, but

1 λ α -32p dCTP

3 λ 3 H dATP

Save \equiv 2' (Same as 2)

For pool spot: Just prior to adding the TdT with the 2 λ of rAn and spot. Then to restore original concentration add 2 λ of H_2O . Then add TdT. (make sure everything well mixed)



3 pool spots

P1

P2 (\equiv Save)

P2

40 minutes

The above isn't such a good way to do this because there is no background level produced from the neg. contr. (I also forgot to save 'save' & acid ppt it by mistake)

1001 T=000.50 A=228682.0(0.7%) B=090116.0(1.0%) C=398834.0(0.5%) R=0.433
 1002 T=000.50 A=010042.0(3.0%) B=650738.0(0.5%) C=792653.0(0.5%) R=5.040
 1003 T=000.50 A=005706.0(5.0%) B=442938.0(0.5%) C=506555.0(0.5%) R=8.071
 1004 T=000.50 A=230364.0(0.7%) B=287778.0(0.7%) C=675748.0(0.5%) R=1.249
 1005 T=000.50 A=015394.0(3.0%) B=842970.0(0.3%) C=141420.0(0.3%) R=5.447
 1006 T=000.50 A=009122.0(3.0%) B=681974.0(0.5%) C=854122.0(0.3%) R=5.777

A

$$\frac{dNTP}{3'} = \frac{.06mM}{4.44 \times 10^5 \mu M} = 1351$$

$$2 \Rightarrow \frac{10012}{(1534)(12.5)} (1351) = 71$$

$$2' \Rightarrow \frac{5706}{(9122)(12.5)} (1351) = 69$$

$$C \quad \frac{dNTP}{3'} = \frac{.20}{4.44 \times 10^5} = 4,505$$

$$2 \Rightarrow \frac{650738}{(842870)(12.5)} (4,505) = 278$$

$$2' \Rightarrow \frac{442938}{(681974)(12.5)} (4,505) = 224$$

$$B \quad \frac{dNTP}{3'} = \frac{.26}{4.44 \times 10^5} = 5,85$$

$$\frac{228682}{230364(12.5)} (5,856) = 465$$

$$I \quad \frac{dNTP}{3'} = \frac{.164}{4.44 \times 10^5} = 0.449$$

$$\frac{99116}{287778(12.5)} (14,414) = 397$$

	#	g	expected
A	70	5.90%	5
C	251	21.2%	25
G	468	39.5%	25
T	397	33.5%	45
	1186	100.1%	

A = .059 C = .211 G = .395 T = .335

term = 1.678049E-02
 ala = .083345
 arg = 9.392546E-02
 asn = 1.900626E-03
 asp = 1.272453E-02
 cys = 7.224946E-02
 gln = 5.651846E-03
 glu = 1.058047E-02
 gly = .156025
 his = 6.797154E-03
 ile = 1.195783E-02
 leu = .1216352
 lys = 1.580374E-03
 met = 7.807176E-03
 phe = 6.127485E-02
 pro = .044521
 ser = 8.340953E-02
 thr = .012449
 trp = 5.226838E-02
 tyr = 1.079169E-02
 val = .132325

d ATP
 B
 C
 A
 G

TdT Tailing

Goal: In the large scale tailing it would be more convenient to use more concentrated DNA & dNTP, while maintaining the same dNTP/3'OH ratio. This has the advantage of a smaller total rxn volume and the requirement for less TdT. Therefore, this TdT assay will use greater concentrations.

Also, it would be nice to compare tailed vs. untailed plasmid on an agarose gel. Some of the DNA from this rxn will be run on a gel.

Procedure

1. 2 λ PUC13/11^{III} @ 1 μ g/1 λ

~~0.5 λ PUC13/11^{III} @ 4 λ per~~

1 λ of 1:5 dil of TdT

2.3 λ 1:40 dil of 50 mM dATP

2 λ 1:10 dil of 50 mM dGTP

6.4 λ 1:10 dil of 50 mM dTTP

2.6 λ 1:10 dil of 50 mM dCTP

1 λ α -³²P dTTP

3.7 λ ³H dGTP

~~0.5 λ PUC13/11^{III} @ 4 λ per~~

25 λ

2 μ g DNA

std

(This makes for 1/10 the amt used in prior assay)

115 mM dATP

4 mM dGTP

1.28 mM dTTP

0.52 mM dCTP

2. Same as above but

1 λ α -³²P dCTP

3.7 λ ³H dATP

25 λ total

1:40 pm

Before adding the TdT, withdraw 2 μ l and spot on disk to count. (Do not replace with 2 μ l H₂O - if the mixture is homogeneous before withdrawal then the concentrations will remain the same; the total volume will just decline slightly.)

After the rxn goes for its designated length withdraw 8 μ l from each rxn mixture ($\frac{8}{25} \times 1 = .64 \mu\text{g DNA}$) to run out on gel. Then stop rxn and count & ppt as before.

Remember, in the calculations, only $\frac{23-8}{25} (2 \mu\text{g}) = 1.2 \mu\text{g}$ will remain

1	226	T=000.50	A=012040.0(2.0%)	R=011006.0(3.0%)	C=023354.0(2.0%)	R=0.356
2	227	T=000.50	A=017036.0(2.0%)	R=025300.0(1.5%)	C=074220.0(0.5%)	R=2.481
P1	228	T=000.50	A=122006.0(1.0%)	R=526060.0(0.5%)	C=261003.0(0.3%)	P=3.795
2	229	T=000.50	A=222076.0(0.7%)	R=370006.0(0.2%)	C=534942.0(0.3%)	S=1.795

Remained only 237-82 = 155

026

was
counted

A

$$\frac{dNTP}{3'OH} = \frac{.115 \mu M}{8.88 \times 10^{-5} \mu M} = 1295$$

$$\frac{18086}{293,886} (1295) = 11$$

C

$$\frac{dNTP}{3'OH} = \frac{.4 \mu M}{8.88 \times 10^{-5} \mu M} = 4505$$

$$\frac{45002}{829,946} (4505) = 31$$

T

$$\frac{dNTP}{3'OH} = \frac{1.28 \mu M}{8.88 \times 10^{-5} \mu M} = 14,414$$

$$\frac{11096}{526,960} (14,414) = 40$$

G

$$\frac{dNTP}{3'OH} = \frac{.52 \mu M}{8.88 \times 10^{-5} \mu M} = 5,856$$

$$\frac{12940}{138,996} (5,856) = 73$$

90

A 11 7%

C 31 20%

G 73 47%

T 40 26%

100%

Remained

1 2 Puc 13/4

A = .07 C = .2 G = .47 T = .26

term = .018382
ala = 9.399999E-02
arg = .111766
asn = .002254
asp = .015134
cys = .056212
gln = 7.560001E-03
glu = .017766
gly = .2209
his = .00644
ile = 9.645999E-03
leu = .088504
lys = .002646
met = 8.553999E-03
phe = .031096
pro = .04
ser = .067134
thr = .014
trp = .057434
tyr = .008372
val = .1222

Goal: - find out min. volume of enzyme
 - get DNA for gels
 - get DNA for HPLC

1. 17 PUC13/HIII @ 1 μ g/17

47 Pout-

1.27 1:40 dil of 50 mM dATP
 17 1:10 dil of 50 mM dCTP
 3.27 1:10 dil of 50 mM dTTP
 1.37 1:10 dil of 50 mM dGTP
 17 α -32P dTTP
 37 α -H dGTP
 -TAT
 -H₂O
 257 total

1 μ g
 STD

.06 mM dATP
 .20 mM dCTP
 .64 mM dTTP
 .26 mM dGTP

~~Q12 TAT undil~~
~~Q32 H₂O~~

Rxn went for 1 1/2 hours

2. same as 1 but,

17 α -32P dCTP
 37 α -H dATP

(A)

17 TdT undil
 8.37 H₂O

(B)

57 1:10 dil TdT
 4.37 H₂O

(C)

37 1:10 dil TdT
 6.37 H₂O

(D)

27 1:10 dil
 7.37 H₂O


Withdraw 27 from the rxn mix before adding
 TdT and spot.

~~After the rxn ends withdraw 137 and save for~~
~~HPLC~~

labelled: A-1 B-1 C-1 D-1
 A-2 B-2 C-2 D-2

A total of 127 was withdrawn from ^{each} A-1 & A-2 (67
 on gel, 67 saved for HPLC) Nothing was withdrawn
 from B, C, or D

A-1 A-2 B-1



I gaged: B-2 & A-2 were washed through the filter before they were acid ppt'd \Rightarrow very few counts were absorbed

1. A = H-dwip, b = 1d1
2. A = 34-dATP, B = 422pdcip

028

1076 T=000.50 A=00000.0(2.0%) B=155591.2(1.0%) C=153972.4(1.07%) D=3.208
1077 T=000.50 A=015000.4(3.0%) B=022404.8(2.0%) C=056230.0(1.57%) D=2.276
1078 T=000.50 A=000100.0(2.0%) B=000000.0(1.5%) C=113114.3(0.87%) D=2.712
1079 T=000.50 A=000000.0(1.5%) B=100000.6(0.77%) C=317160.0(0.77%) D=2.592
1080 T=000.50 A=012724.9(3.0%) B=040464.7(1.5%) C=070590.0(1.57%) D=2.592
1081 T=000.50 A=100000.0(1.0%) B=000000.0(1.0%) C=042673.6(0.77%) D=0.581
1082 T=000.50 A=000000.0(3.0%) B=000000.0(2.0%) C=047170.0(1.57%) D=0.581
1083 T=000.50 A=000000.4(1.5%) B=150000.0(1.0%) C=071778.5(0.77%) D=0.581
1084 T=000.50 A=000000.0(1.5%) B=000000.0(0.5%) C=071778.5(0.77%) D=0.581
1085 T=000.50 A=100000.0(0.77%) B=000000.0(0.5%) C=071778.5(0.77%) D=0.581
1086 T=000.50 A=000000.0(1.5%) B=000000.0(0.5%) C=071778.5(0.77%) D=0.581
1087 T=000.50 A=000000.0(0.77%) B=000000.0(0.5%) C=071778.5(0.77%) D=0.581
1088 T=000.50 A=000000.0(1.5%) B=000000.0(0.5%) C=071778.5(0.77%) D=0.581
1089 T=000.50 A=000000.0(0.77%) B=000000.0(0.5%) C=071778.5(0.77%) D=0.581
1090 T=000.50 A=000000.0(1.5%) B=000000.0(0.77%) C=071778.5(0.77%) D=0.581
1091 T=000.50 A=000000.0(0.77%) B=000000.0(0.5%) C=071778.5(0.77%) D=0.581

Rxn A

(A) $\frac{(16902)(2)}{(92833)(13)} (1351) = 18$ ✓ unreliable due to gas

(C) $\frac{(38495)(2)}{(646835)(13)} (4505) = 42$ ✓ unreliable due to gas

(E) $\frac{(27340)(2)}{(41480)(13)} (5856) = 591$

(I) $\frac{(106581)(2)}{(329711)(13)} (14414) = 717$

Rm B

$$\textcircled{A} \frac{66,58^+}{\cancel{35,160}(2)} (1351) = 3 \checkmark \text{unreliable good}$$

$$\frac{\cancel{40,000}(23)}{24,030}$$

$$\textcircled{C} \frac{1192,833(2)}{(813,210)(23)} (4,505) = 93 \checkmark \text{unreliable good}$$

$$\textcircled{E} \frac{(25,160)(2)}{(49,750)(23)} (5856) = 258$$

$$\textcircled{F} \frac{(69,675)(2)}{(406,142)(23)} (14,414) = 215$$

Rm C

$$\textcircled{A} \frac{15,156(2)}{274,033(23)} (4,505) = 79$$

$$\textcircled{E} \frac{107,580(2)}{214,000(23)} (1351) = 59$$

$$\textcircled{G} \frac{4,734(2)}{22,560(23)} (5856) = 224$$

$$\textcircled{I} \frac{10,464(2)}{26,280(23)} (14,414) = 190$$

Rxn D

$$\textcircled{A} \frac{72,084 (2)}{316,971 (22)} (1351) = 27$$

$$\textcircled{C} \frac{(158,686) (2)}{(711,080) (23)} (4505) = 87$$

$$\textcircled{G} \frac{(14,430) (2)}{(41,504) (23)} (5856) = 177$$

$$\textcircled{T} \frac{(25,521) (2)}{(197,240) (23)} (19,414) = 162$$

Results	Mg enzyme		Slap enzyme		Bla enzyme		Zn enzyme	
	Rxn A		Rxn B		Rxn C		Rxn D	
dNTP	#	q ₀	#	q ₀	#	q ₀	#	q ₀
A	18*		3*		59	10.7	27	6.0
C	42*		93*		79	14.3	87	19.2
G	591		258		224	40.6	177	39.0
T	717		215		190	34.9	162	35.8
		sum: 1308		sum: 473	552	100.6	453	100.6

* unreliable due to goof up

Razor
slit
orientation
slit

dNTP? →

6 hour autoradiogram of the above gel

PvuII HindIII Digest

1907 PvuII at 100/17 (= 1907)
 33 ~~487~~ HindIII at 40/17 (= ~~487~~ 0)
 26.47 10X "core" buffer 142
 264.47 total

157
 H₂O

Ran went O/N @ 37°C

Given 200µg from Ray - 10µg in 107) Swed

19/10/07
 old cut
 new cut
 A-1

cut by
 A-1

After 12 hours the digest was stopped. Together w/
A-Z (from the day before) the DNA was ~~ethanol~~ EtOH.

The ~~PCR~~^{Pac13/HTII} was re-eluted in 140 μ TE

The A-Z was re-eluted & cut w/ ~~EcoRI~~ ~~BamHI~~ ~~SmaI~~
EcoRI
at 20 μ / μ

1 μ EcoRI at 20 μ / μ (= 20 μ)

~~EtOH~~ ppt & dry A-Z DNA

2 μ 10X High Salt buffer

1 μ H₂O

20 μ total

checking DNA concentration

	<u>13/H III - batch 1</u>	<u>13/H III - batch 2</u>
280 →	.29	.138
260 →	.378	.277

diluted 1:40 (57 DNA + 1957 TE)

~~40x~~ $\frac{40x}{22} = \frac{\text{mg}}{\text{ml}}$ ^{260 absorbance}

$(\frac{260}{280})$ ratio should be about 2, which it is)

batch-1 .7 $\mu\text{g}/\lambda$ } tail length
 batch-2 .5 $\mu\text{g}/\lambda$ } ~1.5x as long as target

To make 1ml total of the above pool - use 5 times the quantities

- ✓ 500 μ l K^+ cacodylate
- ✓ 100 μ l $CoCl_2$
- ✓ 200 μ l 1:100 dil of dilt
- ✓ 60 μ l 1:10 dATP
- ✓ 20 μ l dCTP
- ✓ 64 μ l dTTP
- ✓ 26 μ l dGTP
- ✓ 30 μ l H_2O

1000 μ l (And this will be a 5X pool \oplus)

dATP is
indicated

Add $CoCl_2$ last (at least after dilt)

TdT Assay

- ✓ 1. 5 λ P α ++(5 λ)
- ✓ 1.42 μ l 13/H^{II} - batch 2 at 0.5 μ g/ λ
- ✓ 0.5 λ TdT
- ✓ 1 λ α -³²P dTTP
- ✓ 2 λ ³H dGTP
- ✓ 15.1 λ H₂O
- 25 λ total

2. same as #1, but

- ✓ 1 λ α -³²P dCTP
- ✓ 2 λ ³H dATP

Take 2 λ out to spot before adding TdT

→ 2 λ (taken out of 24.5)

12:05 pm

8 λ withdrawn from both 1 & 2
before being stopped

⇒ 15 λ total counted

$$\text{pvc's m.w.} = (649)(2686 \text{ g/mol}) = 1,743,214 \frac{\text{g}}{\text{mole}}$$

~~1 x 10⁻⁶ g~~

$$1 \times 10^{-6} \text{ g} \left(\frac{1}{1,743,214} \frac{\text{mole}}{\text{g}} \right)$$

$$1 \mu\text{g} = 5.74 \times 10^{-13} \text{ moles}$$

12.3'OH per molecule

$$\Rightarrow 1.15 \times 10^{-12} \text{ 3'OH M.}$$

in 1 μg :

Scint Data on next page

$$\text{A} \quad \frac{dATP}{3'OH} = \frac{.06 \mu M}{(1.7)(4.44 \times 10^{-5} \mu M)} = 1,931$$

$$\frac{19,380(2)}{(24,300)(15)} (1,931) = 110 \quad \frac{63,290(2)}{167,258(15)} (1,931) =$$

$$\text{C} \quad \frac{dCTP}{3'OH} = \frac{.120 \mu M}{(1.7)(4.44 \times 10^{-5} \mu M)} = 6,435$$

$$\frac{241,348(2)}{511,642(15)} (6,435) = 405$$

$$\text{E} \quad \frac{dGTP}{3'OH}$$

$$\frac{.26 \mu M}{(1.7)(4.44 \times 10^{-5} \mu M)} = 8,366$$

$$\frac{(19,380)(2)}{(24,300)(15)} (8,366) = 476$$

$$\text{F} \quad \frac{dTTP}{3'OH} = \frac{.64 \mu M}{(1.7)(4.44 \times 10^{-5} \mu M)} = 20,592$$

$$\frac{39,378(2)}{193,762(15)} (20,592) = 558$$

	#	%	(um) p25 white
A	97	6.3%	5.9%
C	405	26.4%	21.1%

53 T=000.50 A=010380.0(3.0%) B=039373.0(2.0%) C=552.00
54 T=000.50 A=063290.0(1.5%) B=241343.0(0.7%) C=365.00
55 T=000.50 A=024300.0(2.0%) B=193762.0(0.7%) C=282380.00
56 T=000.50 A=167258.0(0.7%) B=511642.0(0.5%) C=3703.00

Carry
over

{ - counts/prole
{ - at most 200 counts/prole



1
2
13/4
13

A = .063 C = .264 G = .31 T = .363

term = 1.561953E-02
ala = .08184
arg = 8.912469E-02
asn = 2.488563E-03
asp = 1.224531E-02
cys = 7.055631E-02
gln = 6.203736E-03
glu = 7.28469E-03
gly = .0961
his = 1.042826E-02
ile = 1.577961E-02
leu = .1449818
lys = 1.480437E-03
met = 7.08939E-03
phe = 8.261916E-02
pro = 6.969601E-02
ser = .1080773
thr = .016632
trp = .0348843
tyr = 1.433886E-02
val = .11253

razor
index

44 hour
autoradiogram

TdT assay w/ less enzyme / pilot for large scale prep

✓ 1. 5 μ l post^{tt} (5x) std
✓ 1.4 μ l Puc 13/H^{III} - batch 2 at 0.5 μ g/ μ l .7 μ g
1 μ l α -32 P dCTP
2 μ l ³H dGTP
1 μ l 1:3 dil of TdT
✓ 14.1 μ l ~~total~~ H₂O
25 μ l total

2. Same as above, but
1 μ l α -32 P dCTP
2 μ l ³H dATP

- Take 2 μ l out to spot before adding TdT
 \Rightarrow (2 μ l taken out of total volume of 24 μ l)

- Withdraw 8 μ l from each (x2) before stopping
rxn \Rightarrow 15 μ l counted

rxn went for 1/2 hour

001 T=000.50 A=062396.0 (1.5%) B=025949.0 (2.0%) C=104022.0
 002 T=000.50 A=026034.0 (2.0%) B=052517.0 (1.6%) C=099021.0
 003 T=000.50 A=239512.0 (0.7%) B=263518.0 (0.7%) C=556428.0
 004 T=000.50 A=175398.0 (0.7%) B=273325.0 (0.5%) C=584198.0

$$A \quad \frac{26,084 (2)}{175,398 (15)} (1931) = 38$$

$$C \quad \frac{54,510 (2)}{373,320 (15)} (6,435) = 121$$

$$B \quad \frac{62,396 (2)}{239,512 (15)} (8,336) = 290$$

$$T \quad \frac{26948 (2)}{268,618 (15)} (20,592) = 275$$

9


A 38 5.2%

C 121 16.7%

B 290 40.0%

T 275 38.0%

Rec 13
Rec 14
Rec 15



~~The unit is should be ~3x as~~
~~concentrated, but it does not look it. Ask~~
~3x as much cut added as unit

13/01/11
21



Large Scale Random TdT Tailing

175/1.4 = 125x scale

175 ~~175~~ Puc13/H^{III} at 0.5 μ g/l
625 λ 5X pool⁺⁺
50 λ TdT ($\frac{1}{25}$ of the 17 in 252)

87.5 μ g DNA
S¹id

3125
to
1562.5
ecc

~~225 λ 5X pool⁺⁺~~
~~175 λ Puc13/H^{III}~~
~~175 λ Puc13/H^{III}~~

1252 Puc13/H^{III} at 0.5 μ g/l (80)

✓ 87.5 λ Puc13/H^{III} at 0.5 μ g/l

43.75

✓ 312.5 λ 5X pool⁺⁺

S¹id

✓ 25 λ TdT ($\frac{1}{25}$ of the 17 in 252)

✓ 10 λ α -³²P dATP

✓ 20 λ ³H dGTP

✓ 1107.5 λ H₂O

1562.5 λ total

2. Same as above, but

10 λ α -³²P dATP

20 λ ³H dATP

- Before adding TdT withdraw ~~10 λ~~ 10 λ from each tube and spot

- After rxn ends (1 hour) withdraw 20 λ from each tube - place ~~in~~ in culture tube - for counting

- withdraw 7 λ and place in eppendorf - this is for gel

- The remaining bulk should be chilled / placed extracted / ~~extract~~ pp^t.

Did I use the true concentration
of the dNTPs in this calculation
or did I just take them to be 500M

$T=000.50$ $A=015440.0(3.0\%)$ $B=000.00$
 $T=000.50$ $A=011106.0(3.0\%)$ $B=000.00$
 $T=000.50$ $A=153372.0(1.0\%)$ $B=000.00$
 $T=000.50$ $A=000000.0(0.0\%)$ $B=000.00$

$T=000.50$ $A=140868.0(1.0\%)$ $B=360858.0(0.5\%)$ $C=420.00$


$$\boxed{A} \quad \frac{11,106}{140,868(2)} (1931) = 76$$

$$\boxed{C} \quad \frac{32,122}{360,858(2)} (16,435) = 286$$

$$\boxed{G} \quad \frac{15,440}{153,372(2)} (8,336) = 420$$

$$\boxed{I} \quad \frac{5,142}{141,982(2)} (20,592) = 373$$

	#	%
A	76	6.6%
C	286	24.8%
G	420	36.4%
T	373	
total	1155	32.3%
		100.1%



1 2 3 4

A = .065 C = .248 G = .364 T = .323

term = 1.664903E-02
ala = 9.027199E-02
arg = .1004221
asn = 2.412475E-03
asp = 1.350986E-02
cys = 6.713361E-02
gln = 6.91548E-03
glu = 1.015014E-02
gly = .132496
his = 9.20452E-03
ile = 1.335282E-02
leu = .1248611
lys = 1.812525E-03
met = 7.64218E-03
phe = 5.957186E-02
pro = .061504
ser = 9.361386E-02
thr = .01612
trp = 4.279621E-02
tyr = 1.198815E-02
val = .117572

#1 eluted m ~ 110A
#2 eluted m ~ 80A

	Sample	
λ (nm)	1	2
260	.104	.269

$\leftarrow 1/100$ dilution (27 in 198 λ TE)

280 .054 .160

320

$$\frac{40\%}{22} = \frac{mg}{\mu l}$$

$\swarrow A_{260}$

$$1 \Rightarrow \frac{100}{40(.104)} (100) = 18.9 \frac{mg}{\mu l} = \frac{mg}{\mu l}$$

\swarrow dil. factor

$$2 \Rightarrow \frac{100}{40(.269)} (100) = 49 \frac{mg}{\mu l} = \frac{mg}{\mu l}$$

Since these solutions are visibly darker
 \Rightarrow re-precipitate

dilute to total volume of 200 μ l for both 1 & 2

After 2nd phenol extraction / ~~200~~ ppt

	1	2
260	.60	.49
280	.28	.23
320	0	.04

1:200 dil

$\lambda = 260$ reading

$$\boxed{1} \Rightarrow \frac{200(.60)}{22} = 5.5 \frac{\text{mg}}{\text{ml}} = \frac{\mu\text{g}}{\lambda}$$


(165 μg total
in 300)

$$\boxed{2} \Rightarrow \frac{200(.49)}{22} = 4.5 \frac{\text{mg}}{\text{ml}} = \frac{\mu\text{g}}{\lambda}$$

(135 μg total
in 300)

above in a total volume of ~ 300 per
1 μg per 2

1 2 13/4 III



Following 2 ~~0000~~ phenol/ETOH
treatments (should be same
as on p. 37 white, but swirly is not!)

Puc19 check of concentration

260	.353
280	.183
320	0

dilution 1:100

$$\frac{100(.353)}{22} = 1.6 \mu\text{g}/\lambda \quad (\sim 65\lambda \text{ total})$$

Digestion of 19 w/ HindIII

Set aside 5 λ uncut in a separate vial
in original DNA binary vial (96 μg)

60 λ DNA @ 1.6 $\mu\text{g}/\lambda$ (= 96 μg)	
50 λ HindIII @ 10U/ λ (= 500 units)	2.54 μg of PBR rec'd
12.2 λ 10X core buffer	1.52 μg
12.2 λ	

incubated at 37°C

@ 10:10am

After 4 $\frac{1}{2}$ hours an aliquot (0.5 λ) was run
out on a gel

cut uncut

Puc 19
19/III

after
4 hrs \Rightarrow



PXII allowed to continue for 10 hours
then phenol/ETOH'd and recultured in
total of 600 TE (some spilled)

13-uncut
14-uncut
19/III-

Puc 19/III

260

.118

280

.053

320


0

1:200 dil

$100 (.118) = .54 \mu g/l$

Rec 13
Rec 14
Rec 15

040



~~The object is should be ~3x as~~
~~concentrated, but it does not look it. Ask~~
~3x as much cut added as meat

TdT Random Tailing Pilot Rn w/ Pvc 19 (similar to p. 34 yellow)

- ✓ 5 μ l pool⁺⁺ (5x)
- ✓ 1.3 μ l Pvc 19/H^{III} at .54 μ g/1 μ l
- [1 μ l ³²P dCTP
- 2 μ l ³H dATP
- 1 μ l 1:2.5 dil of TdT
- ✓ 14.7 μ l H₂O
- 25 μ l total

std
7 μ g

2. Same as above but

- 1 μ l ³²P dCTP
- 2 μ l ³H dATP

- Before Adding TdT Spot 1 μ l of the rxn mixture
- After Rn is over withdraw 5 μ l and place in separate culture tubes - these will be counted - the remaining 19 μ l should be chilled (& stored for purification methods)


10:45 am \rightarrow 11:30 am
 \Rightarrow 45 minutes

Unfortunately, I forgot to change the filter between 1 & 2, so these counts were combined onto a single filter. About the only thing I can conclude is that the rxn worked ~~OK~~

001° T=000.62 A=056764.5(1.5%) B=099917.7(1.0%) E=191355.4(2.0%) 760
 001° T=000.20 A=052200.0(2.0%) B=065355.0(2.0%) C=143750.0(1.5%) 751
 0220 T=000.20 A=062800.0(2.0%) B=153300.0(1.5%) C=280320.0(1.5%) 751

Of the 197 remaining, 37 each was swept
 (for an unimpacted gel) and the other 167
 was phenol extracted/~~each~~ ppt'd
 and re-eluted in 127 each (per 127)
 DE

1000 2000 1 2 13/11 A2



NZK as much as the DNA sample as
non-extracted sample added (even less
currentated relative to non-extracted than
gel would suggest)

TdT Purification Scheme Trial

Goal: To make TdT to test numerous methods of purification following 18A

- 1. ~~25A~~ pool ++ (5X) std
- ✓ 7A Puc13/H^{III} at 0.5 µg/l 3.5 µg
- 1.5A TdT
- 1A α-³²P dTTP
- ✓ 6.5A H₂O
- 100A T₀4

↑ incubate at 37°C for 1/2 hr

This should have been 125A total ⇒ bases more concentrated (Should not make too much of a diff since this exp. is aimed at separating labeled strand from dNTPs and Encomp. is not to be measured other than visualization on a gel)

25A withdrawn & worked up several BRL prep



After EtOH ppt

reeluted in 10A TE - 6A of this run out on gel

6A Nacs purified 13/III tail (25A → Nacs → 10A)

6A Puc13/H^{III}

1A Puc13/H^{III}

Nas
lawsonii
13/11/44



Much of the "Nas" floated away upon loading because there was still plenty of ~~EtOH~~ not yet dried from the pellets. (And there was a low ratio of loading buffer)

The "Nas" was completely cold upon loading \Rightarrow all the dNTP's were gone (By contrast the phenol/~~EtOH~~ samples are still pretty hot \Rightarrow dNTP's present)

Preparation of NACS

1 gram of NACS
10 ml 2M NaCl TE

- mix gently for 1 hour
- allow mixture to settle for 1 hour - remove supernatant by cap
- repeat until supernatant is clear (none more time)
- store at 4°C until

25% of the previous day's sample was loaded in a 100 syringe Nax column, eluted in 45% ~~0.5M~~ ppt, redissolved in 25%. 67 (25% loading buffer) of the Nax was compared at an equal amount (at equiv. concentration) on a gel

— "Nax"
— "raw"

↳ gel no DNA recovered, at least none appeared on gel (no bands for the picture)

Large Scale Random Tail w/ Pvc 19

42.3% scaled

557 Pvc 19/H^{III} at .54 μ g/l 29.7 μ g
 2127 SX pool⁺⁺ std J
 6.97 ~~222~~ TdI
 - 1.000
 - 4.20
 10587 Total

1. 27.57 Pvc 19/H^{III} at .54 μ g/l 14.9 μ g
 1067 SX pool⁺⁺ std
 8.57 TdI
 57 α -³²P d ATP
 107 γ -H d ATP

✓ 3727 ~~8000~~ H₂O
 10587 total
 5297

2. Same as above but
 57 α -³²P d CTP
 107 γ -H d ATP

1 hour at 37°C

- Before Adding TdI withdraw 57 from each tube & spot as the pool

- After rxn ends (1 hour) withdraw 107 from each tube to save for counting

- The rest is gel loaded for purif.

13/III¹²16/IV¹²

2



193 T=000.50 A=017552.0(3.0%) B=004162.0(5.0%) C=024556.0(2.0%)
194 T=000.50 A=007704.0(5.0%) B=029448.0(2.0%) C=043378.0(1.0%)
195 T=000.50 A=120550.0(1.0%) B=037418.0(1.0%) C=257118.0(0.5%)
196 T=000.50 A=063906.0(1.5%) B=290984.0(0.7%) C=451868.0(0.5%)

[A] $\frac{7704}{63906(2)} (1931) = 116$

[C] $\frac{29448}{290984(2)} (6435) = 326$

[G] $\frac{17552}{120550(2)} (8336) = 607$

corr [T] $\frac{4162}{87418(2)} (20592) = 490$

	#	%
A	116	7.5%
C	326	21.1%
G	607	39.4%
T	490	31.8%
	1539	99.8%

A = .075 C = .212 G = .395 T = .318

term = 2.063025E-02
ala = 8.374001E-02
arg = 9.766374E-02
asn = 2.98125E-03
asp = 1.570125E-02
cys = .0665733
gln = 7.473001E-03
glu = 1.392375E-02
gly = .156025
his = 8.426999E-03
ile = 1.442925E-02
leu = .1149443
lys = 2.64375E-03
met = 9.42075E-03
phe = 5.359572E-02
pro = .044944
ser = 8.311725E-02
thr = .0159
trp = 4.961595E-02
tyr = .0126405
val = .12561

Big gel

← #1

← #2

Purified according to Ray's Protocol, #1 & #2
each eluted in ~130 λ TE

Px13/1100
1:11-tails
2/



<u>Wm</u>	1	2
260	.165	.150
280	.122	.100
320	.036	.017

1:40 dilution

$$1 \Rightarrow \frac{40(.165)}{22} = .3 \mu\text{g}/\lambda$$

~~Ray's Protocol~~
~~5/2/19~~

WJ100000

another way to do calc.
can off people

$$11,552(52.9) = 928,500$$

total
counts put in

$$1833(\overset{90}{100}) = \overset{90}{183,300}$$

144,640

$$\frac{293,290}{928,500} \Rightarrow \frac{1}{3} \text{ of counts received}$$

$$\frac{144,640}{928,500} \Rightarrow \sim 16\%$$

2 λ from both 132 was dotted on filter,
dried & counted

³H³²P

$$T=0.0100, A=0.01533 \cdot 0(5.0\%), B=0.00399 \cdot 0(1.5\%), C=0.02620 \cdot 0(5.0\%), D=0.00000$$

$$T=0.0100, A=0.01533 \cdot 0(7.0\%), B=0.004520 \cdot 0(3.0\%), C=0.007654 \cdot 0(3.0\%), D=0.00000$$

2 λ counted

From counts taken on sample following tailing 1X11

$$\frac{1}{10\lambda} \text{ } ^3\text{H} \Rightarrow \frac{17,552 \text{ cpm}}{10\lambda} = \frac{1,755 \text{ cpm}}{\lambda}$$

$$\text{The above} \Rightarrow \frac{1,833 \text{ cpm}}{2\lambda} = \frac{917 \text{ cpm}}{\lambda}$$

$$\frac{917}{1,755} = 0.52 \times \text{as concentrated}$$

$$\frac{32\text{p}}{10\lambda} \Rightarrow \frac{4162 \text{ cpm}}{10\lambda} = \frac{416 \text{ cpm}}{\lambda}$$

$$\frac{400 \text{ cpm}}{2\lambda} = \frac{200 \text{ cpm}}{\lambda}$$

$$\frac{200}{416} \approx 0.5 \times \text{as concentrated as original sample}$$

$$\text{orig. sample had } \frac{14.9 \mu\text{g}}{52.47} = .028 \frac{\mu\text{g}}{\lambda}$$

$$\text{following gel extraction} = \frac{.028 \mu\text{g}}{2 \lambda} = .014 \frac{\mu\text{g}}{\lambda}$$

(1607 contains 2.24 μg total)

$$\boxed{2} \text{ } ^{32}\text{P}: \frac{7704 \text{ cpm}}{10 \lambda} = 770 \frac{\text{cpm}}{\lambda}$$

$$\frac{1543 \text{ cpm}}{2 \lambda} = 771 \frac{\text{cpm}}{\lambda}$$

$\frac{771}{770} \approx$ the same concentration as original sample

$$\text{ } ^{32}\text{P}: \frac{29,448 \text{ cpm}}{10 \lambda} = 2,944 \frac{\text{cpm}}{\lambda}$$

$$\frac{4529 \text{ cpm}}{2 \lambda} = 2,265 \frac{\text{cpm}}{\lambda}$$

(4) about the same concentration as original PM mix

orig concentration

$$\frac{14.9 \mu\text{g}}{5292} = .028 \mu\text{g}$$

So, purified DNA $\sim .028 \mu\text{g}$
 (1602 \Rightarrow 4.5 μg total)

The DNA was once again phenol
 extracted & EtOH ppt
 - each 1 & 2 was redissolved in 202 TB

17 was dated on filter & counted

New Puc 19

~~26.72 Puc 19 @ 13.07 mg/ml~~

~~1007 ~~26.72~~ Hind III @ 10 u/l~~

26.72 Puc 19 DNA @ 13.07 mg/ml (= 349 μ g)
 1007 ~~26.72~~ Hind III @ 10 u/l (= 1000 units)
 142 10x cut buffer
 140.77 total

37°C @ 9:30 pm

← after
12 hours



Total above rxn mix:

existing { 26.7 Puc19 @ 13.02 mg/ml (349 μ g)
 100.7 HindIII at 10 u/l (1,000 units)
 14.7 10x core buffer
 140.7 total

150.7 H₂O
 50.7 HindIII
 23.9 10x core buffer


new 364.6

total


+180.7 H₂O
 20.7 10x core

544.6 total @ 11 am 7-25-84

Shows
 following
 dilution
 & addition
 of core
 buffer



100
 9400



152 of Lypk. (reconst.) HIII was added
 to Pan Mix $(152 - \frac{94}{2}) = 1354$

Placed into 2 tubes

tube 1

N2407

~~2000~~

(muzzugant)

tube 2

~~2000~~ ~ 2407

(muzzugant)

tube 3

~~1750~~

(muzzugant)

~~For both pt (1) 0.04 1.04/1.04 3.04/1.04~~
~~(extracted by 2.00 assumed equal volume)~~

The mixture was $\phi\phi\phi$, $\phi\phi\phi/\text{CHCl}_3$, and CHCl_3 extracted but not $\phi\phi\phi$ ppt'd. The organic phase was ~~boiled~~ passed and "back extracted"

nm ^{19/11/11}
~~260~~ 500
~~280~~ 278
~~320~~ 001

"back extracted"

305

197

009

diluted 1:40

11

40x

22

$$\text{Concentration } \frac{40(5)}{22} = \text{~~1.8~~}$$

$$= \boxed{91 \mu\text{g/l}}$$

$$(\sim 200 \mu\text{g} \Rightarrow 152 \mu\text{g})$$

$$\frac{260}{280} \text{ ratio } 1.8$$

$$\frac{40(305)}{22} = \boxed{55 \mu\text{g/l}}$$

$$(\sim 220 \mu\text{g} \Rightarrow 121 \mu\text{g})$$

$$1.5$$

TdT Assay

Goal: To check activity of P-L TdT
(based on assay on p. 41-white)

✓ 57 pool⁺⁺ (5X) std
✓ C.87 DNA Pol I/HindIII @ .9/μg/λ 72 μg
2.57 P-L TdT
17 ³²P dCTP
✓ 27 ³H dATP
✓ 13.77 H₂O
257 total

Same as above except
17 ³²P dCTP
27 ³H dATP


45 minute rxn

~~Ex 900~~

~~Ex 900~~

Before adding TdT 27 from each is spotted
After rxn 67 from each is withdrawn & placed
on a gel

051



2-1-54

→

$T=000.50$ $A=073740.0(1.5\%)$ $B=015170.0(3.0\%)$ $C=102.7$
 $T=000.50$ $A=017342.0(3.0\%)$ $B=054534.0(1.5\%)$ $C=090$ $136.0(1.0\%)$ $R=0.205$
 $T=000.50$ $A=199944.0(0.7\%)$ $B=050312.0(1.5\%)$ $C=31$ $252.0(1.0\%)$ $R=3.144$
 $T=000.50$ $A=092355.0(1.0\%)$ $B=257644.0(0.7\%)$ $C=45$ $3612.0(0.7\%)$ $R=0.251$
 $6130.0(0.5\%)$ $S=0.25$

[4]

$$\frac{17,342(2)}{92,356(17)} (1931) = 43$$

[C]

$$\frac{54,534(2)}{257,644(17)} (6,435) = 160$$

[6]

$$\frac{\cancel{199944} 73,740(2)}{199,944(17)} (8,336) = 362$$

[T]

$$\frac{15,170(2)}{50,312(17)} (20,592) = 730$$

	#	%	
A	43	3.3%	$\Rightarrow 1.87$
C	160	12.4%	
G	362	28.0%	
T	730	57.4%	

.quit A = .033 G = .123 G = .28 T = .564
** ED RUN DEASE II **

term = 1.103602E-02
ala = .03444
arg = 3.733212E-02
asn = 7.43143E-04
asp = 6.34703E-03
cys = .1084911
gln = 1.270467E-03
glu = 2.09212E-03
gly = .0784
his = 2.786533E-03
ile = 1.340064E-02
leu = .1689361
lys = 3.40857E-04
met = 5.211341E-03
phe = .218532
pro = .015129
ser = 7.571900E-02
thr = .004059
trp = .0442171
tyr = 1.270644E-01
val = .15792

A new batch of pool⁺⁺ was prepared according to recipe on p. 32w.

TdT Assay: To test out new pool in anticipation of yet another large scale tailing rxn

① 5 μ l pool ⁺⁺ (5X)	Std
✓ 0.8 μ l Puc 19/HindIII @ 91 μ g/l	72 μ g
1.8 μ l P-L TdT	
1 μ l α - ³² P dATP	
2 μ l ³ H dATP	
✓ 13.2 μ l H ₂ O	
24.3-25 μ l total	

2) as above except
1 μ l α -³²P dATP
2 μ l ³H dATP

rxn went 45 min.

Before adding dT spot 2 μ l

After rxn on window for gel

run out on a gel - looked as expected;
not photographed

8. T=000.50 A=002346.0(2.5%) B=001862.0(1.0%) C=0048900.0(2.0%)
 9. T=000.50 A=0011614.0(3.0%) B=018044.0(1.0%) C=030044.0(2.0%)
 10. T=000.50 A=317204.0(6.7%) B=0046506.0(1.5%) C=322394.0(0.5%)
 11. T=000.50 A=152126.0(1.8%) B=11308.0(0.77%) C=450404.0(6.5%)

A) $\frac{11,614 (2)}{21,204 (17)} (1931) = 17$

C) $\frac{13,044 (2)}{211,308 (17)} (6,435) = 47$

G) $\frac{22,346 (2)}{217,204 (17)} (8,336) = 101$

T) $\frac{1,868 (2)}{46,506 (17)} (20,592) = 97.3$

	#	%
A	17	9.1028
C	47	25.39
G	101	54.36
T	21	11.39
	186	100%

	#	%
A	17	6.5%
C	47	17.9%
G	101	38.5%
T	97.3	37.1%
	262	100%

TdT Assay

This assay uses the old "batch 2" of TdT (see p. 3W)

1. 5 λ pool⁺⁺ (5 λ)
- 0.8 λ Pucl9/H^{III} @ 91 μ g/ λ
- 1.5 λ "batch 2" TdT
- 1 λ α -32 P dATP
- 2 λ γ -32 H dGTP
- 14.2 λ H₂O
- 25 λ total


std
72 μ g

2. Same as above except

1 λ α -32 P dATP
2 λ γ -32 H dATP

- Before adding TdT spot 2 λ
- After 1 hr concludes with row 6 for gel check

10 pm 4.5 in 14.



149 T=000.50 A=022346.0(2.0%) B=001868.0(2.0%) C=028900.0(2.0%) S=0.033
 152 T=000.50 A=011614.0(3.0%) B=013044.0(1.0%) C=030044.0(2.0%) S=0.033
 153 T=000.50 A=217204.0(0.7%) B=046506.0(1.5%) C=322394.0(0.5%) S=0.213
 157 T=000.50 A=152126.0(1.0%) B=211303.0(0.7%) C=450404.0(0.5%) S=0.213

A) $\frac{11,614 (2)}{18,204 (17)} (1931) = 17$

C) $\frac{13,044 (2)}{211,308 (17)} (6,435) = 47$

G) $\frac{22,346 (2)}{217,204 (17)} (8,336) = 101$

T) $\frac{1,868 (2)}{46,506 (17)} (20,592) = 97.3$

	#	%
A	17	9.10%
C	47	25.39%
G	101	54.36%
T	21	11.39%
	186	100%

	#	%
A	17	6.5%
C	47	17.9%
G	101	38.5%
T	97.3	37.1%
	262	100%

TdT Assay: To find the proper amt of batch 2

- 1.5 λ pool⁺⁺ (5X-batch 2)
- 0.8 λ pUC19/HindIII @ .9 μ g/ λ
- 6.5 λ "batch 2" TdT
- 1 λ α -³²P dTTP
- 2 λ ³H dGTP
- 9.7 λ H₂O
- 25 λ total

std
.72 μ g

2. Same as above except

- ✓ 1 λ α -³²P dCTP
- 2 λ ³H dATP

- Before adding TdT spot 2 λ

- After 100 minutes ~~antidote~~ for gel check

~~100 minutes~~

100 minutes

001 T=000.50 A=008006.0(3.0%) B=001485.0(1.0%) C=010522.0(3.0%) R=0.185
 002 T=000.50 A=009356.0(3.0%) B=012138.0(3.0%) C=015086.0(3.0%) R=1.297
 003 T=000.50 A=120616.0(1.0%) B=042258.0(1.5%) C=178078.0(0.7%) R=0.365
 004 T=000.50 A=174738.0(0.7%) B=284736.0(0.7%) C=320492.0(0.5%) R=0.365

$$\boxed{A} \quad \frac{9356}{174,738 (11.5)} (1931) = 9$$

$$\boxed{C} \quad \frac{12138}{284,746 (11.5)} (6435) = 24$$

$$\boxed{G} \quad \frac{8006}{120,616 (11.5)} (8336) = 48$$

$$\boxed{T} \quad \frac{1486}{44,258 (11.5)} (20592) = 60$$

	#	%
A	9	6.49%
C	24	17.0%
G	48	34.0%
T	60	42.6%
	141	100%

There appeared to be a lot of glycerol in rxn mixture - maybe so much so that the rxn was inhibited

TdT Assay

1. ✓ 5 λ pool⁺⁺ (batch 2)
- ✓ 0.8 λ Puc19/H^{III} @ .91 μ g/ λ .728 μ g DN_r
- 3 λ "batch 2" TdT
- 1 λ α -³²P dTTP
- 2 λ ³H dGTP
- ✓ ~~13.2~~ 13.2 λ H₂O
- 25 λ total \rightarrow actually 23 λ total since 2 λ withdrawn

2. Same as above except

1 λ α -³²P dCTP
2 λ ³H dATP

45'

Before adding TdT withdrawn 2 λ

\rightarrow withdrawn for gel

\rightarrow checked, not photographed - taken
gear - no detectable exo activity

→ 0.1 NO calculation for full length is wrong. It depends on the concentration of dNTP that is changing in the original assay since path. However, more DNA was added. Revised # of dNTPs/5' out are below.
(See p. 324 back)

51 ng/μl

~~51.051 (1.15 × 10⁻¹² moles per μl)~~

$$59 \mu\text{g} \Rightarrow 55 (1.15 \times 10^{-12} \frac{\text{moles}}{\mu\text{g}})$$

$$= 6.325 \times 10^{-11} \text{ moles of 3' ends}$$

$$\frac{6.325 \times 10^{-11} \text{ moles of 3' ends}}{.0010725 \text{ L}}$$

$$= 5.9 \times 10^{-8} \text{ M}$$

$$.06 \text{ mM} = .06 \times 10^{-3} \text{ M} \\ = 6 \times 10^{-5} \text{ M}$$

$$\frac{\#A \quad 6 \times 10^{-5} \text{ M} \quad \text{dATP}}{\#C \quad 5.9 \times 10^{-8} \text{ M} \quad 3' \text{ ends}} = 1,017$$

$$\frac{2 \times 10^{-4} \text{ M}}{5.9 \times 10^{-8}} = 3,390$$

TdT random Tailing - large Scale

27837 total

121 fold scan

121 λ Puc 19/H^{III} @ .91 μ g/ λ
~~632 λ Puc 19/H^{III} @ .91 μ g/ λ~~
~~632 λ Bacteriophage~~

40 μ g DNA

121 λ Puc 19/H^{III} @ .91 μ g/ λ 110 μ g DNA

429 total scale up

~~429 total~~
~~632 λ Puc 19/H^{III} @ .91 μ g/ λ~~
~~632 λ Bacteriophage~~

1. ✓ 60.5 λ Puc 19/H^{III} @ .91 μ g/ λ 55 μ g DNA
 ✓ 215 λ Puc 19/H^{III}
 129 λ TdT - "batch 2" (@ 100 μ g/ml)
 6 λ α -32 P dCTP
 10 λ 3 H dATP
 ✓ 652 λ H₂O
 1072.5 total

2. Same as above except
 6 λ α -32 P dCTP
 10 λ 3 H dATP

from each tube

- Before adding TdT withdraw 5 λ and spot
- At conclusion withdraw from each tube
 - 10 λ for counting - in culture tube
 - 6 λ for minigel

40 minutes

Note - water was withdrawn
before the stable quantity of
TdT was added \Rightarrow

$$\frac{1072.5 - 2.5}{1072.5} =$$

$$\frac{1072.5}{1072.5 - 2.5} = 1.14 \times \text{as concentrated}$$

$$\text{water } 57(1.14) = 5.687$$

of pool
equiv. drawn

$T=000.50$ $A=003318.0$ (5.0%) $B=000202.0$ (2.0%) $C=004100.0$ (5.0%) $S=0.99$
 $T=000.50$ $A=001248.0$ (1.0%) $B=002002.0$ (2.0%) $C=007135.0$ (7.0%) $S=1.25$
 $T=000.50$ $A=032962.0$ (1.0%) $B=033252.0$ (2.0%) $C=011308.0$ (1.0%) $S=0.45$
 $T=000.50$ $A=036228.0$ (2.0%) $B=076976.0$ (1.5%) $C=074308.0$ (1.5%) $S=0.85$

$$\boxed{A} \quad \frac{1248 \text{ (5.68)}}{36228 \text{ (10)}} (1931) = 38 \quad .0195$$

$$\boxed{C} \quad \frac{2008 \text{ (5.68)}}{76976 \text{ (10)}} (6,435) = 95 \quad .0148$$

$$\boxed{G} \quad \frac{3318 \text{ (5.68)}}{82962 \text{ (10)}} (8,336) = 189 \quad .0227$$

$$\boxed{T} \quad \frac{302 \text{ (5.68)}}{23352 \text{ (10)}} (20,592) = 151 \quad .00735$$

A	38	8.0%
C	95	20.0%
G	189	40.0%
T	151	32.0%
	<u>473</u>	<u>100</u>

249

$\frac{1-87}{P(TAA) (0.8)(.32)(.32) = .00819}$
 $P(TAG) = (.32)(.08)(.40)$
 $P(TGA) = .01024$
 $= .08432$

$$(98.1568) \xrightarrow{100} 15.6\%$$

1914 12



I'm somewhat skeptical of the counting
results - the gel shows a definite
Stubble insert

~~500~~ ~~ug~~ withdrawn & phos extracted
 (1 phenol, 12.11 phenol / chloroform, ether)

#1 phenol'd & checked on a gel-looked good
 So #2 also phenol'd, but not checked

both #1 & #2 need to let ether further
 evaporate from them

~ 1/2 oz #1 loaded into Naes

(~ 27.5 μ g DNA)

→ 520 λ / diluted to 5.2 ml w/ Indigo buffer

eluted in ~~750 μ l~~ 3 x 250 μ l ~~Naes~~

825 λ 3 x 275

short ppt

pellet size recovered quite small

Attempted purification scheme: EtOH ppt -
but first chelate w/ EDTA

100 λ of rxn #1	5 μ g DNA in 1 mM CoCl_2
557 μ of 500 mM EDTA	26 mM EDTA
117 3 M NaCl ^{Sodium} acetate	283 mM NH_4 acetate
2917 95% EtOH	

* Note: Not NH_4 !

pelleted & resuspended in 157 TE

original
- 1000
- 1000

OD measurements

	raw rxn #1 (phenol'd)	EtOH'd #1
260	.439	.738
280	.279	.513
320	0	.057

1:40 dilution

Counting Experiment

37 of #1 raw and 37 of #1 EtOH'd were
dropped on filters dried & counted
"PI" "PIE"

37 of #1 raw and 37 of #1 EtOH'd were
acid ppt and filtered & counted according
to the normal polymerase assay

1 refers to rxn mixture #1

37
acid ppt

A: 34

B: 32p

$$095 \quad T=000.50 \quad A=000586.0(5.8\%) \quad B=010928.0(3.0\%) \quad C=058588.0(7.5\%) \quad P=0.845$$

$$096 \quad T=000.50 \quad A=00175.0(1.0\%) \quad B=021394.0(2.0\%) \quad C=174514.0(3.7\%) \quad P=0.140$$

$$097 \quad T=000.50 \quad A=001090.0(1.0\%) \quad B=000066.0(0.1\%) \quad C=001424.0(1.5\%) \quad P=0.055$$

$$098 \quad T=000.50 \quad A=004544.0(5.0\%) \quad B=000322.0(0.3\%) \quad C=014704.0(3.0\%) \quad S=0.035$$

acid ppt: 37

(these vertical lines
are here by mistake)

Interpretation & calculations from the above

I. First of all, these results (PI & I) should be in agreement with the earlier pool & acid ppt counts

A. P.58W \Rightarrow

$$PI: \left\{ \begin{array}{l} \frac{3}{5}(82,962) \approx 44,566 \\ 34 \quad 49,777 \approx 44,566 \end{array} \right. \Rightarrow \text{counting of pool} \quad \text{for rxn #1 good}$$

$$32p \left\{ \begin{array}{l} \frac{3}{5}(23,352) \approx 10,928 \\ 14,011 \approx 10,928 \end{array} \right. \nearrow$$

B.

$$I: 34 \left\{ \begin{array}{l} .3(3318) \approx 1090 \\ 995 \approx 1090 \end{array} \right. \Rightarrow \text{counting of acid ppt counts for rxn #1 good}$$

$$32p \left\{ \begin{array}{l} .3(302) \approx 66 \\ 10.6 \approx 66 \end{array} \right.$$

III. Recovery on EtOH ppt

(note the phenol extraction is assumed to have 100% - this should be true)

For "1" 3H of the crude mixture had 1090^{3H}
for "1E" 3H had 4644 3H counts

but 1E represents 100% of "1" resuspended
in only 15% TE
//

$$\frac{100}{15} = 6.7 \text{ fold concentration}$$

$$\frac{4644}{1090(6.7)} = 64\% \text{ recovery}$$

Based on 32P counts

$$\frac{328}{66(6.7)} = 74\% \text{ recovery}$$

⇒ ~70% recovery of DNA

IV. 6.9 ^αNTP ~~3H~~ removed by EtOH ppt:

$$3H \left\{ \begin{array}{l} \text{counts due to } \alpha\text{NTP alone} \\ \rightarrow \cancel{44,566(6.7) - 1090(6.7)} \leftarrow \text{before EtOH ppt} \\ \rightarrow \cancel{133,176 - 4644} \leftarrow \text{after EtOH ppt} \end{array} \right.$$

$$3H \left\{ \begin{array}{l} \text{counts due to } \alpha\text{NTP alone} \\ \rightarrow 133,176 - 4644 \leftarrow \text{after EtOH ppt} \\ \rightarrow \frac{128,532}{6.7} = 19,199 \end{array} \right. \quad \begin{array}{l} \text{all seen} \\ \text{due to } \alpha\text{NTP} \\ \text{remaining} \end{array}$$

~~32P~~

~~21394-328~~

$$^{32}\text{P} \left\{ \frac{21394-328}{(10,928-66)(67)} = .29 \right.$$

$^{3}\text{H} \Rightarrow$ 56% of unincorporated dNTP removed

$^{32}\text{P} \Rightarrow$ 71% of unincorporated dNTP removed

Conclusions

1. ~70% of the DNA was recovered by EtOH ppt. (concentration of DNA computable from this fact)
2. ~60% of the unincorporated dNTP was removed
3. A relatively small percentage of the counts even following EtOH ppt are due to the DNA. Despite removal of over half of the unincorporated dNTP only ~1.5% of the total counts are acid pptable
4. When EDTA was added & the ~~reservoir~~ salt switched to Na acetate (in place of ammonium acetate) degradation of the DNA did not occur
5. Perhaps the Vites purified stuff—although very cold—is just pure DNA. The pellet was not very big, but counting it should tell for sure—after all EtOH ppt of

Ethanol ppt of the remaining raw rxn mix

#1 1072.57 - 6202 = 4537

#2 1072.57

#1

4537 g #1 raw rxn (231 mg of DNA)
 257 g 500 mM EDTA
 507 3M Na acetate

3507 g #1 100% EtOH 18ug of DNA

19.37 g 500 mM EDTA

397 g 3M Na acetate

5177 g EtOH

12257 total

~~507 g #1~~

~~3.37 g 500 mM EDTA~~

~~10.57 3M Na acetate~~

~~2507 g EtOH~~

~~5257~~

#2 - without DNA for gel reference & quantifying
 divide into 5 aliquots of 3507 and treat
 as #1

#1 Phenol:chloroform TE

#2 Phenol:chloroform in 50% TE

2 aliquots for quantifying

code 1EA ← ethanol ppt as opposed to "Peri
 = 100" ← "Sp. 4"

2- control
2- ethanol ppt
1- ethanol ppt -
large bubble



This gel shows that the DNA has maintained
its integrity following EtOH ppt.

27g each counted

4131 cpm remaining

4131 cpm remaining

27g each counted

I. Do these results give earlier pool & acid precipitable count data?

A. p. 58W \Rightarrow check of "pool" results

$${}^3\text{H} \left\{ \begin{array}{l} \frac{2}{5} (36,228) \approx 10,080 \\ 14,491 \approx 10,080 \end{array} \right.$$

\Rightarrow fairly consistent agreement

$${}^{32}\text{P} \left\{ \begin{array}{l} \frac{2}{5} (76,976) \approx 20,348 \\ 30790 \approx 20,348 \end{array} \right.$$

The ${}^3\text{H}$ data should be more accurate for this comparison because the ${}^{32}\text{P}$ is in continuous decay.

B. Check of acid pptable count results

$${}^3\text{H} \left\{ \begin{array}{l} .2 (1248) \approx 254 \\ 250 \approx 254 \end{array} \right.$$

\Rightarrow fairly consistent agreement

$${}^{32}\text{P} \left\{ \begin{array}{l} .2 (2008) \approx 254 \\ 402 \approx 254 \end{array} \right.$$

II. Recovery of EtOH ppt

A. Run #1

$$3H \left\{ \begin{array}{l} \cancel{1030(15)} \\ 3052(15) \\ 1090(17.5) \end{array} \right. \xrightarrow{27} \begin{array}{l} \cancel{1030(15)} \\ 246 \text{ recovery} \end{array}$$

1EA 3507 represents 3507 of 1A

resuspended in only 20x $\Rightarrow \frac{350}{20} = 17.5x$ concentration

$$\begin{array}{l} \text{upon} \\ \text{recounting} \end{array} \xrightarrow{21} 32p \left\{ \begin{array}{l} \cancel{2578(15)} \\ 66(17.5) \end{array} \right. = \begin{array}{l} \cancel{279} \text{ recovery} \\ 3356 \text{ recovery} \end{array}$$

$$.25(18\mu g)$$

$$20x \searrow$$

$$.23\mu g$$

$$20x \rightarrow 4.5\mu g$$

Clearly, something
here is wrong
at 3H & 32p don't give

B. Run #2

$$3H \left\{ \begin{array}{l} 4131 \rightarrow \text{upon recounting} \\ 6794 \\ 254(20) \end{array} \right. = \begin{array}{l} \cancel{1346} \text{ recovery} \\ 816 \text{ recovery} \end{array}$$

#2 was ~10007 redissolved in 507 \Rightarrow 20 fold dilution

$$32p \left\{ \begin{array}{l} 2526 \\ 254(20) \end{array} \right. = 50\% \text{ recovery}$$

\Rightarrow on the average ~ ~~65%~~ recovery
65% recovery

$$\Rightarrow \frac{.65(55\mu g)}{507} = .72\mu g/7$$

$$507 = 36\mu g$$

$$0.01\mu g$$

III. removal of dNTP

A. rxn #2

$$\frac{10,080 - 254}{(37,992 - 6794)}$$

$$34 \left\{ \frac{37,992 - 6794}{(10,080 - 254)20} = .16 \right.$$

⇒ 84% of
all unincorp
dNTP's
removed

$$32 \left\{ \frac{147,462 - 2,526}{(20,348 - 254)(20)} = .36 \right.$$

⇒ 64% of
all unincorp
dNTP's
removed

∴ ~75% of all
unincorp. dNTP's removed

B. rxn #1

$$34 \left\{ \frac{10910 - 66}{(103,540 - 3052)(1.5)} = .19 \right.$$

$$(44,566 - 1090)(17.5) \Rightarrow 80\%$$

all unincorp
dNTP's removed

$$32 \left\{ \frac{(17,252 - 2578)(1.5)}{(10,928 - 66)(17.5)} = .12 \right.$$

⇒ 88% of
all unincorp
dNTP's removed

Attempted Annealing, Fill & ligation

Anneal

2.8 μ l #2 rxn/phenol'd / EtOH'd @ .72 μ g/l

.8 μ l 1M tris pH 7.5

36.4 μ l H₂O

40 μ l total

2 μ g DNA
20mM Tris

1. heat to 68°C for 5 minutes

2. incubate at 58°C for 1½ hours

Fill-in & ligation

✓ 40 μ l of #2 rxn/p'd/e'd DNA at 50 μ g/l (from above) 33.3 μ g/l
 ✓ 12 μ l of 5' γ 14mer PstI from BRL lypk (=60 units) 10 μ l
 ✓ 1.2 μ l 15mM dITP
 ✓ 0.6 μ l 1M ATP
 ✓ 0.6 μ l 1M MgCl₂
 3 μ l of 2U/2 T4 DNA ligase (=6 units)
 ✓ 0.5 μ l 5mM dATP
 ✓ 0.5 μ l 5mM dCTP
 ✓ 0.5 μ l 5mM dGTP
 ✓ 0.5 μ l 5mM dTTP
 59.4 μ l total

in 20mM tris pH 7.5
(=2 μ g)
BRL lypk
10 μ l
42 μ M
42 μ M
42 μ M
42 μ M

3:20pm Aug. 15

incubated at room temp = 22.2°C

at 11am (=6½ hours)


157 colonies - 100% transformed

✓ 500 μ g = 5 μ g

Reaction
annealing
mix

parent
#2-12

derived
mix-5



Is it circular or linear - try
restricting it! yes, but restriction
sites will randomly appear also assuming
all went well

unique sites → NATR
in PstI or BclI → ECKRV
multiplexing 1161
regions

~~1000~~ ~~amplified~~ mix (= 267 ng DNA)

8:45 - 1 ml of overnight \rightarrow 200 ml XGT

The following plates were made

PUC19 (unmodified; mutant)

The PUC19/racchar

0.5 ng

0.5 ng

(from
"amplified
mix")

5 ng

5 ng

50 ng

50 ng

500 ng

The DNA was diluted accordingly

\rightarrow incubated at 2:15 pm

Transformation Protocol

Ca⁺⁺ Cells

1. Pick a single colony (or from stock glycerol) of JM83 and inoculate into a 200ml 2xYT at 37°C in a Shaker overnight. (Use a large flask ~~200ml~~ for good aeration = 300ml)
2. Take 1ml of this overnight culture and inoculate a 200ml 2xYT culture (a large flask again) for 2hr. While growing place a 250ml serum bottle at -20°C to cool.
3. At A₆₀₀ = .6 or "just a little less" harvest culture by centrifugation at 5K for 8 minutes at 0-4°C.
4. Decant supernatant (clear and amber). Resuspend cell pellet in ice-cold 50mM CaCl₂ (filter sterile). Use 1/2 original grown volume (~100ml); resuspend by swirling; do not pipette to resuspend. Incubate on ice for 20 min. w/ occasional swirling.
5. Repellet for 5K at 0-4°C for 8 min.
6. Decant supernatant (clear and colorless) and resuspend cells in 1/10 original volume (~20ml) of ice cold 50mM CaCl₂.
7. Store cells on ice in cold room.

NB:

(i) keep everything cold! (on ice)

For amp selection: add
1ml of 2xYT to each
tube and incubate in 37°C
water bath for 1 hour

transfection

1. For every desired plate, place 0.3ml Ca^{++} cells in a small culture tube and add DNA. Incubate 40 minutes on ice (20 minutes should be sufficient)

Meanwhile, make top agar, 2ml per plate ~~XXXX~~
(2x4T heated in microwave with .9g agar; .09g/l can
don't let boil over in microwave; prepare in small
conical culture flask) and keep warm on hot
plate set at "4.2 and low" (=

2. After the 40 minutes is up heat shock the tubes for 2 min in a $42^{\circ}C$ water bath
(make this bath in a styrofoam tub w/ sink water,
warm to room temp one plate to be used)

3. Remove from heat and add ~~XXXXXX~~ to the soft agar the following on
a per tube amount

per tube, but add to soft agar mix!

15 μ IPTG (1M)

75 μ XGAL (=6, 10 diMeF₅₀)

6 μ 3X amp @ -5 μ g/ml \Rightarrow 75 μ g/ml

4. Add 2ml of the above mix per each tube of Ca^{++} cells; quickly vortex and pour on to 3 agar AMP plates; swirl to coat the whole plate

5. Incubate upside down at 37° overnight

Results:

More or less a lawn on all plates
~~at control~~ even neg. control \rightarrow All blue

A few blue colonies on the periphery of
 the "19" plates "A-50" had 2 or 3 blue
 colonies on the periphery

Making more "anneal II":

- 2x the protocol on p. 65 was followed \rightarrow incubation at \approx to
- also "anneal control" at 1x the original
 50 μ m protocol was made. It was heated only -
 to see if the slower migration is merely
 an artifact of chipping under heat - NO
 Kleenex or ligase added - set out at room temp
 w/ "anneal II"

Exn allowed to go O/N starting 7:15 pm
 \rightarrow 8:30 am
 \downarrow


Bham 1X1

Exon 1

Exon 2

Exon 3

Exon 4



Transformation Data

<u>λB</u>	<u>Stock</u>	<u>[plate]</u>
AMP	25 μ g/ml	50 μ g/ml
tet	12.5 μ g/ml	12.5 μ g/ml
Strap	10 μ g/ml	25 μ g/ml
Chlor	34 μ g/ml	10 μ g/ml

→ The 2^o λ B were not added to the top agar - one in the master plate

puc19

Anneal Mix

DNA
Dilutions

1 μ g/plate
⇒ 3.87 μ g / 1:100,000

5 μ g/plate ⇒
3.87 μ g / 1:10,000

50 μ g/plate
⇒ 38.57 μ g / 1:10,000

[Stock] = 13 μ g/ λ

1 μ g/plate ⇒ 1.57 μ g / 1:100

5 μ g/plate
⇒ 15.7 μ g / 1:100

50 μ g/plate ⇒ 1.57 μ g / 1:100

500 μ g/plate ⇒ 15.7 μ g / 1:100

Media

B broth:

10g bacto tryptone

8g NaCl

add 1 liter dH₂O & autoclave (121°C)

B plates: Add 15g/ λ Agar for plates

Nomenclature for Replica Plating

W \Rightarrow white colony picked from PUC RAN plates

WB \Rightarrow blue colony " " " "

B \Rightarrow blue colony picked from PUC 19 plates

* The number following the above designation ~~refers~~ is arbitrarily assigned to the colony #



W*2 \Rightarrow These colonies were picked from the plate according to the above nomenclature (minus the *), but they were NOT replica-plated. Because I worried for the plasmid prep and ^{plasmid} any cell that would not grow I picked a couple of colonies that would not "be replaced off the loop" by replica plating and were used to inoculate liquid culture (for the plasmid prep) only.

Rice Plasmid prep from Transfected cells 070

1 2 3 4 5 6 7 8 9



The DNA was
eluted in
20x

4x was
run on the
gel

281

281

105

[Handwritten signature]

Smile to face

1 ml 5M NaOH
250 μ l 20% SDS
3.75 ml H_2O
5 ml

found

132g

4 ml 10%

30 ml 129

for 15 and

Hand
1070

135

1.29 NaOH

7.5ml 20% SDS

$$\frac{142.5}{190 \text{ ml S}}$$

200 ml Sol'n III

• 16 g NaOH

1 ml 20% SDS

19 ml 100

20 mil total

.209

1.25 ml 20% ~~SDS~~

23.75 ml

25 ml

$\begin{array}{r} \swarrow 109\text{ lbs} \\ 80\text{ ml} + 240\text{ ml} = 320\text{ ml} \\ 141.92\text{ g} \cdot 64\text{ g} \\ 8\text{ mls} \\ 72\text{ mls} \\ \text{H}_2\text{O} \\ \hline 40 \end{array}$

Soln I = 50 mM Glucose
 25 mM Tris pH 8
 10 mM EDTA (pH 7.5)
 12 mg/ml lysozyme (fresh) (for cosmids make 12 mg/ml). *See over*

Soln II = 0.2 N NaOH
 1% SDS *(keep 1 week) make fresh!!!*
 1.25 ml 10% SDS
 11.25 ml H₂O
 0.19 N NaOH = 1 pellet

Soln III = 5 M KOAC pH 4.8 (60 ml 5M KOAC + 11.5 ml HOAC → 100 ml).

Small total
 200 µl 10% SDS
 100 µl 10N NaOH
 4.65 ml H₂O

*1 ml NaOH
 5% SDS
 H₂O*

*Small
 2 ml water
 1 ml 10% SDS
 2 ml H₂O*

- Spin down 1 ml of o/n culture in Epp tube 15 secs.
 Remove all medium w/ 1 ml Gilson. (If treating large nos. of samples, can use with drawn out pasteur pipette on a water-pump. Ditto for rest of protocol).
See over for 5 ml. cosmid cultures).

*→ 50% w/o lysozyme
 then add 50% w/ 2x lysozyme*

- Suspend pellet in 100 µl Soln I with Eppendorf & vortexing.
 (for cosmids, transfer to Eppendorf tube)
 5 min at R. Temp.
- Add 200 µl Soln II with Eppendorf. Mix by rocking tube sharply but don't vortex.
 Inc. on ice for 5 min.
- Add 150 µl precooled Soln III. Seal cap well, invert tube, mix gently on lid, vortex gently briefly, invert tube (to normal orientation); vortex gently briefly.
- Inc. on ice/5 min.
- Spin 1 min/microfuge.
- Transfer supernatant (360 µl) to fresh tube, avoiding sediment.
- ~~Add 720 µl EtOH~~ *DEP → 65° for 10-15 min?*
- Leave 2 min at R. temp.

- Spin 1 min/microfuge.
- Remove supernatant, carefully but completely with 1 ml Gilson.
- Wash with 500 µl 70% EtOH, rinsing sides of tube with inclining tube (no vortexing or spinning)
 Remove carefully but completely with 1 ml Gilson.

OK

Dry (2 min/high vacuum).
 Take up in 50 µl TE, suspending pellet w/ Eppendorf.
 Vortex, 4' RT, vortex.

*GENOL EXTRACT
 MECP*

Before removing aliquots for digestion (~2 µl) (5 µl for cosmids)
 Clarify by spinning 10 sec/microfuge.

42
 Add RNase in digests. *Can be scaled up preparatively after adding KOAC, cooling & spinning.*
 Use 0.5 M NaOH.

Restriction Digests of Quick plasmid prep

37 of plasmid prep DNA
1.57 μ g Tag I restriction Endonuclease @ 10 u/ μ l
17 μ g 10x core buffer
4.57 μ g H_2O

incubate at 37°C


+ 19 stock digest


.37 μ l Puc 19
1.57 μ g Tag I
17 μ g 10x core buffer
7.27 μ g H_2O

40 min at 37°C

Taq 7 digests

W1-W8-W9-W5-B1-B2-W61-W62-19






cut

control

Electrophoresis results




cut


control

34 32 33 34 35 36 37 38 39 40 61 62 P418 2/10/85


UNCUT



CUT




073



Uncut

Standards from 7/Hind III
9/2/19/Tag I



<u>W1</u>	<u>W2</u>	<u>W3</u>	<u>W4</u>	<u>W5</u>	<u>W6</u>	<u>W7</u>
1.1 Kb	1.3 Kb	1.3	1.3	Ø	1.9	2.2 Kb
2.5 Kb	.85	1	.8		.6	1.3 Kb
2.0 Kb			.6			.5 Kb
1.3 Kb						
1.1 Kb						

(short & fast calculation
of size of bands)

growth after 15 hours in culture, where [X] tried:

fox-5 1.98/ml 1.97 stock

fox-6 8.58/ml .97 stock

	<u>gulf</u>	<u>trimeter</u>	
ST-5	5.8/ml	1.8/ml	357 g 1:10 stc
ST-6	2.58/ml	.58/ml	1.87 g 1:10 stc

FOX4-A has growth (but FOX3-A₃ does not, while all other foxes do) except for FOX4-19 & FOX3-19. (From 9-25 → 9-27 there were allowed to sit at room temp.)

3 (1 ML) ^{paras} of FOX4-A was withdrawn & pelleted; SU decanted & pellet frozen at -70°C (Stored in Ray's "In Use" box) to save for future plasmid prep.

100 ml of FOX4-A was dumped onto a 2x4T plate and α B sensitivity discs applied.

From the 2^o plates (STX & FOX) all colonies were lifted (one at a time w/ a loop) and added to 2.5 ml liquid 2x4T, stored overnight at 4°C, and then 100 ml of each was dumped onto a 2x4T plate

Selection Experiment

Goal: To see if any of the PUCRAN DNA will confer drug resistance

Methods:

JM83 was transformed with PUCRAN DNA = "AII" @ 300 ng per plate (total of 2 plates) according to the protocol on p. 66. (Control: pUC19 @ 5 ng) (AII \rightarrow ~75% white colonies; pUC19 ~100% blue)

"lifts" \rightarrow On

6 ml of 2xYT liquid was added to plate - swished around and raked (w/ glass rake). 200 μ l of this infective liquid was plated onto large diam. petri AMP, 2xYT plates (AMP @ 50 μ g/ml). Antibiotic disks were added to the plate (~30 different ones). There were a total of 4 plates: 2 "lifts" from pUC19 controls (in Pakove) and 1 each from the two AII (= PUCRAN) DNA. (The II merely denotes the second batch prepared according to the protocol on p. 65).

diameter of de ~~1.8~~ ~~1.7~~ ~~1.6~~ ~~1.5~~ ~~1.4~~ ~~1.3~~ ~~1.2~~ ~~1.1~~ ~~1.0~~ ~~0.9~~ ~~0.8~~ ~~0.7~~ ~~0.6~~ ~~0.5~~ ~~0.4~~ ~~0.3~~ ~~0.2~~ ~~0.1~~

		<u>Pu19#1</u>	<u>Pu19#2</u>	<u>ATI#1</u>	<u>ATI#</u>
0.0000					
gofranboin	F/M	2cm	2.2cm	2cm	1.9
olfamethazol	TH	Ø	Ø	Ø	2
lfamethazole	trimethoprim SXT	1.1cm	1.6cm	Ø	
mpicillin	TM	Ø	Ø	Ø	1
tracycline	TE	2cm	2.3cm	2.1	2
polymixin B	PB	1.2cm	1.2cm	1.2	1.2
	NN	1.1cm	1.1cm	1.1	1.1
cefuroxime	UT	Ø	Ø	Ø	Ø
	CF	Ø	Ø	Ø	Ø
ticarcillin	TIC	Ø	Ø	Ø	Ø
gentamicin	GM	1.2cm	1.3	1.1	1.2
moxalactam	MOX	2.2cm	2.7	2.9	2.5
	FOX	2.4cm	2.4	2.4	2.4/1.2cm
cefotaxime	CTX	2.5cm	2.6	2.7	2.6
piperacillin	pip	Ø	Ø	Ø	Ø
0.0000 cefepime	CEP	Ø	Ø	Ø	Ø
	A.V	1.5cm	1.5	<u>1.4</u>	1.4
chloramphenicol	C	2.7cm	2.7cm	2.4	2.4/1.2

Sulfamethoxazole + trimethoprim = SXT
 cefoxitin = fox 30

ST-A colonies picked
from single
colonies growing
around SXT discs
on 2° plates

ST-19 colonies } picked
* EX-19 colonies }
from (and)
on big 19-plates

EX-1 colonies picked
from single
colonies on 2°
plates grown
in absence of
fox (but
incubated
first)

Liquid Growth

Cefoxitin stock

$$50 \text{ mg/ml} \\ = .05 \text{ g} \times 10 \text{ ml} = .5 \text{ g in } 10 \text{ ml}$$

Each tube: 5 ml 2x4T
10% of 25 mg/ml AMP \Rightarrow 508/ml

Cefoxitin (AMP is labelled 10 MCB from disc package; FOX is labelled 30 mcg)

fox-1 37.58/ml

fox-2 758/ml

fox-3 1508/ml

fox-4 3008/ml

3.875 g stock
7.5 g stock

157 g stock

307 g stock

~~6000000000~~

Sulfamethoxazole trimethoprim

stock: 16 mg/ml - trimeth
80 mg/ml - sulf

Sulf Trimeth
~~500 mg~~ ~~250 mg~~

ST-1

108/ml

28/ml

.77 g stock

ST-2

208/ml

48/ml

1.37 g stock

ST-3

408/ml

88/ml

2.57 g stock

ST-4

808/ml

168/ml

57 g stock

for λ b disc sensitivity testing.

In going through old plates a single colony was noticed on a month old TET plate, it was picked & grown 2 days ~~in~~ in air shaker at 37°C (in 12.5 μmol , 6.3 μmol , 3.2 μmol , and 1.6 μmol [TET] against controls)

For all puc19 controls, a loop was dipped into an area of the lawn on the very large, 1" λ b sensitivity plates (away from the discs)

ligand
Colonies picked for plasmid preps

TET @ 1.6 μmol + 50 μmol AMP	2 days 37°C	
19 @ 50 μmol AMP	2 days 37°C	2 days RT
+ 50 μmol sulfa	2 days RT	TET trans plate
18 μmol trimeth		

SXT @ 50 μmol AMP	2 days 37°C	selection
+ 50 μmol sulfa	2 day RT	2" plate
18 μmol trimeth		

SXT @ 50 μmol AMP	2 days 37°C	Selection
50 μmol sulfa	2 days RT	20 plates
16 μmol trimeth		

FOX @ 50 μmol AMP	2 days 37°C	uncol
200 μmol cotribi	2 days RT	2 plates

Pursuit of SXT resistance

colony isolation

The SXT^{res} colonies from the original screening were plated onto 2^o amp/km¹² (12 grid plates). All colonies from this 2^o grid plate were picked and inoculated into 2ml 2x liquid. 100 μ l of this inoculum was spun onto a 2xYT plate and checked for AMP & SXT sensitivity w/ the ϕ paper discs. A pattern similar to the original screening was observed. (The photo below is a picture of the second screening plate)



← A loop taken from as close to the paper disc as possible was smeared onto an LAMP plate to isolate single colonies. Single colonies were isolated and grown on an LAMP grid plate in which the agar prickles were made around the edge of an SXT ϕ paper discs. These are then isolated single SXT^{res} amp^{res} colonies (Puck prepared as above - colonies picked from the lawn on the original selective plate - did not grow at all near the SXT discs).

see
p. 744

The above single colonies were used in the following liquid inoculum:

Each tube: 5ml 2XYT
5081ml AMP

	Sulf trimeth		
ST-4	100 8081ml	1681ml	57g stock
ST-3	40	8	2.57g stock
ST-2	20	4	1.37g stock
ST-1	10	2	.77g stock
ST-0	0	0	0

Following the label
an "A" denotes positive ST's and a "19" indicates
puc19 control.

The colonies were inoculated into the 2XYT+AMP
and allowed to grow for 1 hour 15 minutes
before the SXT was added.

Puc19 controls are from single colony.
Puc19's on an AMP plate

Results

	A	19	
ST-4	0	0	
ST-3	+	+	
ST-2	+	+	
ST-1	0	+	
ST-0	+	+	0 = NO + = grow

ST-4 A & ST-0 19 tried for plasmid prep

peps digested!

37g DNA

17g 10X Cn buffer

1.57g Tag I 581ml AMP

27g 100μg/ml (100μg/ml stock) RNase

2.57g H₂O

107 total

65°C 2 hours

10X Cn buffer 2.57g

581ml AMP

100μg/ml (100μg/ml stock) RNase

8:35am 2.57g

19/Tag
ST/Tag
SH1/Tag
SH2/Tag
Tag
SH1
SH2

